



All-in-one ICR Swab

Environmental monitoring
in critical areas

ICR Swab

Ready-to-use and easy-to-handle

Products for environmental monitoring within isolators and clean rooms must meet the highest safety standards. This includes reproducible and secure detection of contaminants. The ICR Swab is designed for qualitative microbial presence/absence tests. It is ideal for monitoring surfaces that are inaccessible to contact plates.



Benefits

- Clean room and isolator use
 - Triple-bagged
 - γ -sterilized at 25-35 kGY
 - Low abrasion Swab material
 - H_2O_2 -impermeable inner bags
 - Hanging device integrated
- Safe and easy handling
 - Reduced contamination risk
- Residue-free sampling
 - Swab pre-moisturized with NaCl solution
 - No growth promoting residue on surfaces
- Tube rack with hanger
 - High quality stainless steel
 - Suitable for autoclaving
 - Persistent against H_2O_2 and disinfectants

Description

Product name	Order No.	Packaging Units
ICR – Swabs (Swab tip moisturized with NaCl (sodium chloride); triple-bagged, gammasterilized at >25 kGy; TSB with lecithin + Tween® 80 in a reservoir)	1.46529.0100	100
Tube rack for ICR – Swabs	1.46530.0001	1



Open the Swab and sample the desired surface.



Insert the Swab back into the tube. Snap the reservoir containing the culture broth.



Squeeze the reservoir until the broth has fully saturated the Swab and incubate at the appropriate temperature. Monitor for turbidity.

To place an order or receive technical assistance

In Europe, please call Customer Service:

France: 0825 045 645

Germany: 069 86798021

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: [MerckMillipore.com/offices](https://www.MerckMillipore.com/offices)

For Technical Service visit: [MerckMillipore.com/techservice](https://www.MerckMillipore.com/techservice)

[MerckMillipore.com](https://www.MerckMillipore.com)



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

