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Solvent Evaporation

In the solvent evaporation technique, polymer and drug are dissolved in an organic solvent with a low boiling point, followed by evaporation and subsequent dehydration. The chosen organic solvent is selective toward one block, which results in the formation of micelles in nonaqueous solutions. The outcome is usually determined by the type of solvent, the concentration of polymer and drug, and the rate of evaporation. The limitation of this approach lies in the limited choice of solvents, and there is no guarantee the resulting particles will be well defined core-shell particles that can be easily re-dissolved in water.

Example Method


1. Dissolve 2 mg of drug and 20 mg of polymer in methanol (or any other low-boiling solvent that can dissolve both components).
2. Evaporate solvent under vacuum.
3. Add distilled water, incubate at 40 °C for 10 min, and vortex to obtain a clear solution.

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Molecules 2015, 20, 20397–20408, Matthew McKenzie et al.

Hydrogel-Based Drug Delivery Systems for Poorly Water-Soluble Drugs

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
Co-solvent Evaporation

Co-solvent evaporation proceeds by adding water directly to the organic solvent to cause the self-assembly of the micelle and encapsulation of the drug. The outcome is controlled by the type of solvent, the ratio between organic solvent and water, the concentration of water and drug, rate of solvent evaporation, and the order and rate of mixing.²⁰ This approach is limited by the choice of solvent, but usually results in higher drug encapsulation efficiencies.

Example Method

1. Dissolve 20 mg of polymer and 2 mg of drug in 1 mL of acetone, THF, or acetonitrile.
2. Add 2 mL of water dropwise to the organic solvent (or vice versa).
3. Mix for 4 h, followed by the evaporation of the organic solvent.

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Dialysis

Dialysis is probably the most versatile and most common technique used for drug encapsulation since it allows the use of high-boiling solvents such as DMSO, which is removed by dialysis and replaced with water. Although this approach is applicable to many solvent systems, the drug loading efficiency is usually lower than the co-solvent evaporation and the technique can be time-consuming.

The slow process can aid the formation of thermodynamically stable morphologies. A final dialysis step to remove solvent and free drug is often crucial to obtain a product free of organic solvent while maintaining maximum drug loading.

However, while extensive dialysis can assist in the thorough purification of the product, it also can cause the release of the already encapsulated drug and result in low drug encapsulation efficiencies.

Example Method


1. Dissolve 20 mg of polymer and 2 mg of drug in 1 mL of DMF.
2. Add 5 mL of water slowly, with the help of a syringe pump, if possible, to control the rate of water addition.
3. Dialyze against water using a tubular cellulose membrane (Sigma Prod. No. Z726176).

参考文献： Journal of Colloid and Interface Science 354 (2011) 116–123, Zhimei Song et al.

Curcumin-loaded PLGA-PEG-PLGA triblock copolymeric micelles: Preparation, pharmacokinetics and distribution in vivo

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Flash Nanoprecipitation

Flash Nanoprecipitation

Flash nanoprecipitation is a relatively new technique that offers a more rapid solution than other time-consuming methods. Fast mixing and precipitation into a non-solvent for the drug and one polymer block results in a kinetically trapped structure. Although the resulting structures do not have well-defined internal phase boundaries, as would be the case in thermodynamically stable structures, the approach provides an alternative to achieve a fast throughput.

Example Method

1. Prepare a solution of 40 mg of polymer and 20 mg of drug in 1 mL of THF.
2. Use a confined impinging jet mixer to mix the solution with 1 mL of water.
3. Introduce the exit stream into 8 mL of water:THF (9:1 v/v%).
4. Dialyze against water using a tubular cellulose membrane.

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本紙記載の製品は試験・研究用です。ヒト、動物への治療、もしくは診断目的として使用しないようご注意ください。なお、品目、製品情報、価格等は予告なく変更される場合がございます。予めご了承ください。記載内容は2020年8月時点の情報です。©2020 Merck KGaA, Darmstadt, Germany. All rights reserved. Merck, the vibrant M, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.