English

Vented Millex®-AA Particulate Filter Unit

Cat. No. SLAAV255F (50/box)

- Single use only
- Sterile
- Non-pyrogenic

Introduction

The Vented Millex®-AA filter unit is a particulate filter for aqueous solutions dispensed with a syringe or dispensing unit. It will remove particles, precipitates, and undissolved powders larger than 0.8 micron (µm). The sterile Vented Millex®-AA filter unit is non-pyrogenic and non-toxic. This single-use product consists of a 0.8 µm membrane filter sealed in a plastic housing. A 0.2 µm hydrophobic filter in the dome prevents air locking by automatically venting air introduced upstream. Typical applications include the clarification of water, buffers, reconstituted antibiotics, ampouled drugs, bacterial toxins, antiseptic and allergenic preparations, and other aqueous solutions where protein binding is not a concern, as well as applications where air venting is required.

Chemical Compatibility


The Vented Millex®-AA filter unit is compatible with most aqueous solutions. It may be used to filter the agents listed below. This guide has been developed from technical publications, materials suppliers, and laboratory tests, and is believed to be reliable. However, because of variability in temperature, concentrations, duration of exposure, and other factors outside of our control which may affect the use of the unit, no warranty is given or it is to be implied with respect to such information. Agents not listed below should be tested with the Vented Millex®-AA filter prior to use.

Acetic acid (5%)	Butyl alcohol	Hypo (photo)	Petroleum based oils
Aliphatic ethers	Freon® solvent, TF or PCA	Isobutyl alcohol	Phenol (0.5%)
Ammonium hydroxide (6 N)	Glycerin (glycerol)	Isopropyl alcohol	Silicone oils
Benzyl alcohol (1%)	Helium	Mineral spirits	
Boric acid	Hexane	Nitrogen	
Brine (sea water)	Hydrogen (gas)	Pentane	

Directions for Use

WARNINGS:

- To ensure sterility, do not use this product if the package is damaged.
- Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter unit and/or personal injury.
- CAUTIONS:
 - Do not resterilize or reuse this filter unit, as Merck Millipore Ltd. cannot assure the sterility, integrity, and performance beyond a single use.
 - Do not use the Vented Millex®-AA unit at temperatures above 45 °C (113 °F).
 - Do not use the Vented Millex®-AA unit to filter emulsions or suspensions.
 - Do not use the Vented Millex®-AA unit to filter 5 mg or less of active drug materials unless binding studies have been performed.

Specifications

Materials	Mixed cellulose ester, type AA, pore size 0.8 µm Hydrophobic PTFE (polytetrafluoroethylene), 0.2 µm PVC (polyvinyl chloride)
Dimensions	Inlet to outlet Diameter Filtration area Temperature limit Pressure limit at 25 °C Filtration volume Hold-up volume Sterilization method Connections Flow rate at 2.1 bar (30 psi), 25 °C
	25 mm (0.98 in.) 29 mm (1.14 in.) 4 cm² (0.62 in²) 45 °C (113 °F) 3.1 bar (45 psi) inlet and differential maximum 1–100 mL ≤ 0.1 mL after air purge Ethylene oxide gas Female Luer-Lok™ inlet, male Luer-slip outlet ≥ 450 mL/min

Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Merck Millipore Ltd. ("Millipore") or an affiliate. Neither Merck Millipore Ltd. nor any of its affiliates assumes responsibility for any errors that may appear in this document.

Technical Assistance

For more information, contact the office nearest you. In the U.S., call 1-800-MILLIPORE (1-800-645-5478). Outside the U.S., go to our web site at www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

Standard Product Warranty

The applicable warranty for the products listed in this publication may be found at www.millipore.com/terms (within the "Terms and Conditions of Sale" applicable to your purchase transaction).

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