

# New Fast and Innovative Detection of Beer Spoilage Organisms

Jvo Siegrist

Microbiology Focus Edition 1.4

## HYBRISCAN® AN INNOVATIVE SCREENING METHOD FOR BEER SPOILAGE ORGANISMS BASED ON THE DETECTION OF rRNA

The popularity of beer remains high but the quality of beer has to be very high to survive in a competitive market. Beer spoilage organisms are either lactic acid bacteria belonging to the genera *Lactobacillus* and *Pediococcus* or they are obligate anaerobes of the species *Pectinatus* and *Megasphaera* (**Table 1**). Within the species of *lactobacilli* known to cause spoilage of beer, only certain strains can grow in the beer and are responsible for spoiling (exception: *Lactobacillus lindneri* all strains cause spoilage). *L. brevis* is the most common beer spoilage bacterium followed by *L. lindneri*.<sup>1</sup> Additionally many wild yeasts are responsible for beer spoilage such as *Saccharomyces cerevisiae* and *Candida pelliculosa*.<sup>2</sup> One of the biggest problems is biofilm formation in beer plants, which makes it very difficult to remove spoilage organisms.

Table 1. Species of beer-spoilage microorganisms that can be detected with the HybriScan®D Beer-Kit (Cat. No. 62533 , 96 assays)

<b>Genus Lactobacillus:</b>	<i>Lactobacillus acidophilus</i>
	<i>Lactobacillus brevis</i>
	<i>Lactobacillus brevisimilis</i>
	<i>Lactobacillus buchneri</i>
	<i>Lactobacillus casei</i>
	<i>Lactobacillus collinoides</i>
	<i>Lactobacillus coryniformis</i>
	<i>Lactobacillus curvatus</i>
	<i>Lactobacillus fermentum</i>
	<i>Lactobacillus fructivorans</i>
	<i>Lactobacillus lindneri</i>
	<i>Lactobacillus malefermentans</i>
	<i>Lactobacillus parabuchneri (frigidus)</i>
	<i>Lactobacillus paracasei</i>
	<i>Lactobacillus paraplantarum</i>
	<i>Lactobacillus plantarum</i>
<i>Lactobacillus rhamnosus</i>	
<b>Genus Pediococcus:</b>	<i>Pediococcus acidilactici</i>
	<i>Pediococcus claussenii</i>
	<i>Pediococcus damnosus</i>
	<i>Pediococcus inopinatus</i>
	<i>Pediococcus parvulus</i>
	<i>Pediococcus pentosaceus</i>
<b>Genus Pectinatus:</b>	<i>Pectinatus cerevisiiphilus</i>



Highly skilled lab staff perform microbiological analysis in specific quality control laboratories. Most of the laboratories still use conventional standard based cultivation methods, which are very time consuming and take 3 to 5 days for beer to be released to the market.

HybriScan® Beer kit, a rapid test system developed for a faster and more reliable product release of beer and could act as an alternative for the detection of beer spoilage contaminants. After as little as two hours (pre-enrichment for 24 h, if necessary) the brewery could have the first reliable results.

A variety of applications have been developed for HybriScan® including the detection of bacteria and yeast in non-alcoholic beverages. The robustness of the HybriScan® assay enables it, in contrast to other rapid test systems, to detect bacterial contamination in brewer's yeast and leads to efficient use of this valuable resource. Furthermore HybriScan® test system is a perfect tool for microbiological control of dispensing equipment. The legal standard for sterility control of dispensing equipment is 100,000 cfu/mL; a fast, direct determination of beer spoiling bacteria is possible without pre-enrichment-procedure delivering results within two hours.

### Comparison of HybriScan® and other rapid test systems

Performing quality control by using the standard cultivation based method takes a long time. In recent years many companies have developed rapid test systems to hasten this procedure. For quality control of beer and beverages three main technologies are available:

- HybriScan® (sandwich hybridization)
- PCR (Polymerase Chain Reaction)
- VIT (Vermicon Identification Technology)

A comparison of these different technologies is given in **Table 2**. Comparing HybriScan® to PCR or VIT-technology the benefits of this rapid test system are:

- Fast and cost efficient analysis
- Inexpensive read-out technology
- High sensitivity and specificity

Using two different probes for detection of microbial RNA, false-positive results are almost impossible. In **Figure 1** results of quantification of *Lactobacillus buchneri* within a starter culture (silage) of three different samples are presented. Comparison of the HybriScan® test with a cultivation based analytical method (MRS agar) displays the equivalent results within the limits of microbiological sample variability.




**Table 2.** Comparison of different technologies for detection of beer spoilage bacteria.

	<b>Cultivation based method</b>	<b>PCR</b>	<b>VIT</b>	<b>HybriScan®</b>
<b>Method</b>	cultivation based method with optical or microscopic read out	PCR/real-timePCR	fluorescence microscopy	sandwich hybridization and photometrical signal read out
<b>Detection spectrum</b>	detection and identification of all beer spoilage microorganisms	identification of all relevant Beer spoilage microorganisms possible	<i>Lactobacillus</i> sp. and <i>Pediococcus damnosus</i> , <i>Lactobacillus</i> sp. + <i>L. brevis</i> , <i>Pectinatus</i> + <i>Megasphaera cerevisiae</i>	identification of all relevant Beer spoilage microorganisms possible
<b>Sample preparation</b>	selective pre-enrichment	enrichment and lysis of bacteria, if necessary pre-enrichment	selective pre-enrichment	enrichment and lysis of bacteria, if necessary pre-enrichment
<b>Time</b>	3 to 7 days	3 hours to 2 days	2 days	3 hours to 2 days
<b>Costs per test</b>	ca. 1 €	12 €	15 €	3 €
<b>Detection limit (cfu)</b>	1	1-5 x 10 <sup>5</sup>	1 x 10 <sup>3</sup>	1-5 x 10 <sup>3</sup>
<b>Devices</b>	None	PCR cycler	fluorescence microscopy	microplate reader
<b>Advantages</b>	high sensitivity, relatively cheap	high sensitivity, quantitative analysis	simple detection technology set up, detects only living cells (RNA)	rapid and sensitive, qualitative and quantitative detection of living cells, cost efficient analysis
<b>Disadvantages</b>	time consuming, no detection of non-culturable microbes, labor expensive	expensive devices needed, no discrimination between live and dead cells, not officially accepted	time consuming, low sample throughput, expensive, not automatable, difficult data analysis, not officially accepted	no differentiation of serotypes or subspecies, limited probe design (rRNA target), not officially accepted

The vibrant M, Millipore, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.  
© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.





**Figure 1.** Lactobacilli (beer spoilage organisms)

## References

1. Storåards E, Suiiko M, Pot B, Vanhonacker K, Janssins D, Broomfield PLE, Banks JG. 1998. DETECTION AND IDENTIFICATION OF LACTOBACILLUS LINDNERI FROM BREWERY ENVIRONMENTS. 104(1):46-54. <https://doi.org/10.1002/j.2050-0416.1998.tb00974.x>
2. Timke M, Wang-Lieu NQ, Altendorf K, Lipski A. 2008. Identity, beer spoiling and biofilm forming potential of yeasts from beer bottling plant associated biofilms. *Antonie van Leeuwenhoek*. 93(1-2):151-161. <https://doi.org/10.1007/s10482-007-9189-8>

