

# Amsphere® A+ Chromatography Resin


A next-generation Protein A affinity resin engineered for superior process efficiency

Protein A capture chromatography efficiently purifies antibodies from complex mixtures of host cell proteins, DNA, media components, and other impurities through highly specific affinity binding interactions.

Amsphere® A+ is a next-generation Protein A resin that offers outstanding dynamic binding capacity and the performance characteristics needed for high titer, intensified processes.

## Key Benefits

- Outstanding dynamic binding capacity optimized for high-titer and intensified processes
- Excellent durability and caustic stability for enhanced bioburden control
- Superior pressure-flow properties and impurity removal ability




## Technical Characteristics

	Amsphere® A+ Resin
Type of Chromatography	Protein A Affinity Chromatography
Base material	Methacrylic copolymer
Average particle size (d50)	~60 µm
Functional group/ Ligand	Recombinant protein A expressed in <i>Escherichia coli</i>
Dynamic Binding Capacity (DBC)	≥60 mg/mL at 4 min residence time (RT) and 10% breakthrough (BT)
Recommended bed height range	5–25 cm
Working pH range	1–13
Shipping buffer and recommended storage buffer	50 mM sodium phosphate buffer containing 16% Ethanol, pH 7.5
Storage temperature	+2 °C to +8 °C

## Outstanding Dynamic Binding Capacity


Amsphere® A+ resin achieves outstanding DBC using an optimized base matrix pore size and controlled ligand conformation and orientation. These enhancements enable the resin to effectively handle the requirements of high harvest titers in intensified processes.



**Figure 1.** DBC of Amsphere® A+ resin and other commercially available resins with three antibody-like molecules.

## Caustic Stability


Improved ligand linking chemistry results in Amsphere® A+ resin having high caustic stability, sustaining over 94% DBC after 24 hours in 1M NaOH. This enables intensified resin cleaning between cycles, enhancing bioburden control.



**Figure 2.** DBC of Amsphere® A+ and other commercially available resins after 24h in 1 M NaOH.

## Pressure vs. Flow Properties

60 µm base beads enable scalable flow rates up to 335 cm/h for large columns at 20 cm bed height and 2 bar pressure drop. Higher flow rates means faster loading which improves productivity for intensified processes.



**Figure 3.** Column pressure against flow rate for different sized columns containing Amsphere® A+ resin.



## Ordering information

Description	Cat. No.
<b>Amsphere® Bulk Resin</b>	
Amsphere® A+ 25 mL	BP-AMS-A+-0025
Amsphere® A+ 50 mL	BP-AMS-A+-0050
Amsphere® A+ 500 mL	BP-AMS-A+-0500
Amsphere® A+ 5 L	BP-AMS-A+-5000
<b>Amsphere® Prepacked Columns</b>	
Amsphere® A+ in MiniChrom® Column 1 mL	BP-AMS-A+-CN01
Amsphere® A+ in MiniChrom® Column 5 mL	BP-AMS-A+-CN05
Amsphere® A+ in RoboColumn® 8 x 200 µL	BP-AMS-A+-CN016
Amsphere® A+ in RoboColumn® 8 x 600 µL	BP-AMS-A+-CN048

For additional information, contact us.



## Chromatography Columns and Systems

Chromatography columns and systems are critical factors to the successful separation of your valuable molecule. From screening to large-scale production, our columns, systems and single-use solutions are designed to provide robust, consistent performance while providing you with the processing flexibility required in today's changing production environment.



To place an order or receive technical assistance:  
[SigmaAldrich.com/support](https://SigmaAldrich.com/support)



For local contact information:  
[SigmaAldrich.com/offices](https://SigmaAldrich.com/offices)

MilliporeSigma  
400 Summit Drive  
Burlington, MA 01803  
[SigmaAldrich.com](https://SigmaAldrich.com)

We have built a unique collection of life science brands with unrivalled experience in supporting your scientific advancements.

**Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®**

© 2026 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M, BioReliance, Millipore, Amsphere, Milli-Q, SAFC, Sigma-Aldrich, and Supelco are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MS\_DS15259EN Ver. 1.0 69370 04/2026

