



RABBIT ANTI-HUMAN PROSURFACTANT PROTEIN B (proSP-B) POLYCLONAL ANTIBODY

CATALOG NO.:	AB3430
LOT NUMBER:	
QUANTITY:	100 µL
SPECIFICITY:	Reacts strongly with human and mouse proSP-B (detects both the N-terminal and C-terminal propeptides). Does not react with mature SP-B. On human tissue immunostains alveolar Type II cells and non-ciliated bronchiolar cells. On mouse tissue immunostains alveolar Type II cells and non-ciliated bronchiolar (Clara cells) in distal and terminal airways.
IMMUNOGEN:	Full length recombinant peptide (381 a.a.) of the human SP-B (pre)proprotein.
APPLICATIONS:	Immunohistochemistry: 1:4,000-1:8,000 on human tissue with and without antigen retrieval (see protocol on back) and 1:1,000-1:2,000 on adult and fetal mouse lung (see notes on back). Immunoblot: 1:1,000-1:5,000 ELISA: 1:500-1:1,000 Optimal dilutions must be determined by the end user.
SPECIES REACTIVITY:	Human and mouse.
FORMAT:	Rabbit serum.
PRESENTATION:	Liquid. Contains no preservative.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to twelve months. Avoid repeated freeze/thaw cycles.
REFERENCES:	Lin, S., et al., <i>J. Biol. Chem.</i> (1996) 271 :19689-19695. Lin, S., et al., <i>Biochim. Biophys. Acta</i> (1996) 1312 :177-185. Morotti, R.A., et al., <i>Human Pathology</i> . (1999) 30 :618-625. Nogee, L.M., et al., <i>Am. J. Respir. Crit. Care Med.</i> (2000) 161 :973-981.

IMMUNOHISTOCHEMISTRY FOR AB3430

IMMUNOHISTOCHEMISTRY

Note: The enzymatic reaction product was enhanced with nickel cobalt to give a black precipitate.

Sections for antigen retrieval were immersed in sodium citrate buffer, pH 6.0, and heated in a microwave for 15 minutes at 90°C. Endogenous peroxidase was quenched for 15 min. at RT with 3% H₂O₂ in methanol. Sections were blocked with 2% normal goat serum in PBS with 0.2% Triton for 2 hours at RT before incubating overnight at 4°C with the properly diluted AB3430. The sections were then washed 6 times, 5 min. each, in PBS with 0.2% Triton and then incubated for 30 min. at RT with biotinylated goat anti-rabbit antibody (such as Chemicon catalog number AP132B) properly diluted in the blocking solution. Sections were incubated with the ABC complex diluted in blocking solution for 30 min. at RT; NiDAB in 0.1M acetate buffer for 4 min. at RT; Tris Cobalt for 4 min. at RT; and counterstained with nuclear fast red for 2 min. at RT. Controls included (1) formalin-fixed, paraffin-embedded, surgical samples of human pediatric and adult lung immunostained at the recommended dilutions or age-matched controls from wild-type, FVB/N, fetal and embryonic mouse pups and (2) omission of the AB3430 to check for endogenous biotin and peroxidase activity, as well as nonspecific binding of the secondary antibody.

Notes:

- 1) Dilutions used are for normal (i.e. uninjured) adult mouse lung and have been used on FVB/N, Balb/c, nu/nu and SCID mice. The use of antigen retrieval systems to enhance immunostaining for this antibody is not necessary in the mouse lung and may actually increase background staining.
- 2) Immunostaining of Clara cells for SP-B and proSP-B in the adult mouse lung can be variable from sample to sample. Most often when immunostaining is present it is seen in the most distal and terminal airways in the lung. It is specific for Clara cells and does not stain ciliated cells. It has not been found in the adult mouse tracheal cells or glands. Balb/c mice and strains of mice derived from a Balb/c background, such as nu/nu mice of SCID mice, appear to have increased amounts of SP-B in the bronchiolar epithelium by immunohistochemistry. (FVB/N = NIH3 < Balb/c.nu/nu < 'American' CB.17 scid/scid < 'European' CB.17 scid/scid mice (Jennings, Am. J. Respir. Cell Mol. Biol. (1995) 13:297-306 and Zsengeller, Hum. Gene Ther. (1997) 8:1331-1344.))

- REFERENCES:** Shi, S-R. et al., *J. Histochem.Cytochem.* (1991) **39**:741-748.
Gown, A.M., *Appl. Immunohistochem.* (1993) **14**:267-274.
Leong, A. et al., *Appl. Immunohistochem.* (1993) **14**:267-274.
Vorbroke, D.K., et al., *Am. J. Physiol.* (1995) 268 (Lung. Cell Mol. Physiol. 13):L647-656
Wert, S.E., et al., *Pediatr. Res.* (1996) **39**:355A.
Wert, S.E., et al., *Pediatr. Res.* (1995) **37**:356A.
Zhou, L., et al., *Dev. Biol.* (1996) **175**:227-238.
Zhou, L., et al., *J. Histochem.Cytochem.* (1996) **44**:1183-1193.
Jain-Vora, S., et al., *Am. J. Respir. Cell Mol. Biol.* (1997) **17**:1-11.
Zsengeller, Z.K., et al., *Hum. Gene Ther.* (1997) **8**:1331-1344.
Morotti, R.A., et al., *Human Pathology.* (1999) **30**:618-625.
Tokieda, K., et al., *Am. J. Respir. Cell and Mol. Biol.* (1999) **21**:463-472.
Nogee, L.M., et al., *Am. J. Respir. Crit. Care Med.* (2000) **161**:973-981.



Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com