Supelco®

www.sigmaaldrich.com

1.14753.0001
1.14753.0007MQuant®
Chloride TestCl⁻

1. Method

Determination with color-disk comparator

Chloride ions react with mercury(II) thiocyanate to form slightly dissociated mercury(II) chloride. The thiocyanate released in the process in turn reacts with iron(III) ions to form red iron(III) thiocyanate. The chloride concentration is measured **semiquantitatively** by visual comparison of the color of the measurement solution with the color fields of a color disk.

2. Measuring range and number of determinations

Measuring range / color-scale graduation	Number of determinations
3 - 6 - 10 - 18 - 30 - 60 - 100 - 180 - 300 mg/l Cl⁻	200

3. Applications

Sample material:

Groundwater, surface water, and seawater (after dilution)
Drinking water and mineral water
Waters from aquaculture
Boiler and boiler feed water, cooling water
Industrial water
Wastewater and percolating water
Process water
Swimming-pool water
Food after appropriate sample pretreatment
Soils after appropriate sample pretreatment

4. Influence of foreign substances

This was checked individually in solutions containing 100 and 0 mg/l Cl⁻. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %					
Ag⁺	1	Cu²⁺	100	NO₂⁻	100
Al³⁺	100	F⁻	100	Pb²⁺	100
Br⁻	1	Fe³⁺	100	PO₄³⁻	100
Ca²⁺	1000	Hg²⁺	10	S²⁻	1¹⁾
Cd²⁺	100	Mg²⁺	100	SiO₃²⁻	1000
CN⁻	1	Mn²⁺	100	Zn²⁺	1000
Cr³⁺	100	NH₄⁺	1000		
Cr₂O₇²⁻	1	Ni²⁺	100		
				NaNO₃	20 %
				Na₂SO₄	0.25 %

¹⁾ In cases of higher concentrations, eliminate sulfide ions by adding hydrogen peroxide (1 drop of Perhydrol® per 10 ml of sample).

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

2 bottles of reagent Cl-1
2 bottles of reagent Cl-2
1 graduated 6-ml plastic syringe
2 test tubes with screw caps
1 color-disk comparator

Other reagents and accessories:

Hydrogen peroxide 30 % H₂O₂ (Perhydrol®) for analysis EMSURE®, Cat. No. 1.07209

MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535

Ammonia solution 25 % for analysis EMSURE®, Cat. No. 1.05432

Nitric acid Titrisol® for 1mol/l, Cat. No. 1.09966

Chloride standard solution Certipur®, 1000 mg/l Cl⁻, Cat. No. 1.19897

MQuant® Flat-bottomed tubes with screw caps for MQuant® tests with color-disk comparator (12 pcs), Cat. No. 1.17988

Refill pack:**Cat. No. 1.18322**

Chloride Test

Refill pack for 1.14753 and 1.14401

(Reagents **without technical accessories** for the number of determinations stated in section 2)

6. Preparation

- Analyze immediately after sampling.
- The pH must be within the range 1 - 12.**
Adjust, if necessary, with ammonia solution or nitric acid.
- Filter strongly turbid samples.

7. Procedure

	Measurement sample right-hand tube (A) behind the color disk	Blank left-hand tube (B) behind the color disk	
Pretreated sample (10 - 30 °C)	6 ml	6 ml	Inject into the test tube with the syringe.
Reagent Cl-1	6 drops ¹⁾	-	Add, close the tube, and mix.
Reagent Cl-2	6 drops ¹⁾	-	Add, close the tube, and mix.

Hold the comparator to the light, keeping it upright, and rotate the disk until the closest possible color match is achieved between the two large windows.

Read off the result in mg/l Cl⁻ shown in the small window.

¹⁾ **Hold the bottle vertically while adding the reagent!**

Notes on the measurement:

- The color of the measurement solution remains stable for at least 60 min after the addition of reagent Cl-2.
 - Turbidity in the measurement solution makes the color comparison more difficult.
 - If the color of the measurement solution is equal to or more intense than the darkest color on the scale, repeat the measurement using **fresh**, diluted samples until a value of less than 300 mg/l Cl⁻ is obtained.
- Concerning the result of the analysis, the dilution must be taken into account:

Result of analysis = measurement value x dilution factor

8. Method control

To check test reagents, measurement device, and handling:

Dilute the chloride standard solution with distilled water to 100 mg/l Cl⁻ and analyze as described in section 7.

Additional notes see under **www.qa-test-kits.com**.

9. Notes

- Reclose the reagent bottles immediately after use.
- Rinse the test tubes and the syringe **with distilled water only**.
- The contents of the test tubes as well as the test reagents must not be run off with the wastewater!**
Information on disposal can be obtained at www.disposal-test-kits.com.

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