


# MultiScreen® Caco-2 Assay System

## protocol note



## ***Optimization of Caco-2 cell growth and differentiation for drug transport assay studies using a 96-well assay system***

Jeanne E. Phillips, Ph.D. and Andrew Arena  
Millipore Corporation, Life Science Division, Danvers, MA


### **Introduction**

Understanding and testing the absorption, distribution, metabolism, excretion and toxicity (ADMET) properties of candidate compounds is critical to successful new drug discovery. The important first step, absorption, is a drug's ability to cross the target cell's membrane barrier from the point of administration to the site of action. This absorption can occur through passive transcellular or paracellular diffusion, active carrier transport or active efflux mechanisms.


Typically, an immortalized cell line is used as a drug absorption model in studies to understand the drug's permeability and absorption mechanisms. A recognized *in vitro* model system for measuring drug absorption uses the Caco-2 cell line, an immortalized, heterogeneous cell line, derived from a human colorectal adenocarcinoma.<sup>1,2</sup>

Millipore has developed a 96-well filter plate system, the MultiScreen Caco-2 assay system (Figure 1), that is automation compatible and allows for the measurement of drug transport in a high throughput format. The

MultiScreen Caco-2 plate system facilitates the use of *in vitro* model cell lines in the measurement of drug absorption rates. Data generated using this device can be used to rank order the oral absorption profiles of new candidate drug compounds. The following protocol provides guidance for the optimization of the growth and differentiation of the Caco-2 cell line on the MultiScreen Caco-2 plate. Representative drug transport results using an optimized protocol are also presented.



**Figure 1:** MultiScreen Caco-2 assay system components (cat. MACA C02 B5, with single-well feeder plate) and 96-well transport analysis plate (cat. MACA C0R S5). Not shown, assay system components with 96-well feeder plate (cat. MACA C02 S5).



**Figure 2:** Drug transport study using MultiScreen Caco-2 assay system and Caco-2 cells.

## Protocol

### I. Medium, Additives and the Cultivation of Caco-2 Cells

The following medium, additives and protocol are recommended for the cultivation of the Caco-2 cell line (cat. HTB-37) purchased from the ATCC® (Manassas, VA). Dulbecco's MEM with high glucose (cat. D5796), non-essential amino acids (NEAA, cat. M7145), HEPES (cat. H0887), penicillin, streptomycin, and L-glutamine (cat. G1146), EDTA (cat. E8008), trypsin/EDTA (cat. T3924) were all purchased from Sigma Chemical Co. (St. Louis, MO). Fetal bovine serum (FBS, cat. 16000044) was purchased from Invitrogen (Carlsbad, CA).

- a. Heat inactivate FBS by incubating for 1 hour in a 50 °C water bath.
- b. Prepare Caco-2 cell culture medium consisting of Dulbecco's Modified Eagle's Medium (DMEM) with 4.5g/L glucose, L-glutamine and sodium bicarbonate supplemented with:
  - 10% heat inactivated Fetal Bovine Serum (FBS)
  - 0.1 mg/mL of streptomycin
  - 100 units of penicillin
  - 10 mM HEPES buffer
  - 1X Non-essential Amino Acids (NEAA)
- c. Cultivate cells in T-75 flasks in a cell culture incubator set at 37 °C, 5% CO<sub>2</sub>, 95% relative humidity. Allow cells to reach 80 – 90% confluence before detaching and splitting. Do not allow cells to become over confluent (>90%), as this will impact subsequent monolayer formation of cells cultivated for drug transport on the MultiScreen Caco-2 assay system.

## Tips & Techniques

### Converting cells/cm<sup>2</sup> to cells/mL or cells/well

#### 24-Well Filter Plates

Surface area: 0.3 cm<sup>2</sup>/well

Volume: 0.4 mL/well

Cells/cm <sup>2</sup>	Cells/mL	Cells/well
82,000	61,500	24,600

#### 96-Well Filter Plates

##### (MultiScreen Caco-2 plate)

Surface area: 0.11 cm<sup>2</sup>/well

Volume: 0.075 mL/well

Cells/cm <sup>2</sup>	Cells/mL	Cells/well
82,000	120,000	9,000

- d. Rinse cultivated cells in T-75 flasks with 5 mL EDTA. Aspirate off, add 1.5 mL trypsin/EDTA, and incubate at 37 °C for approximately 5 to 10 minutes or until the cells detach and float. This can be confirmed by periodic visual inspection of the flasks.
- e. Add 12 mL of Caco-2 cell culture medium to detached cells and pipette up and down to disperse cells. Dispense 2.5 mL into 6 new T-75 flasks. Add 12.5 mL of Caco-2 cell culture medium for a total of 15 mL per flask and replace in cell culture incubator.

## II. Optimization of Caco-2 Cell Seeding Density on MultiScreen Caco-2 Assay System

Though Caco-2 cells commercially available for this application originate from the same immortalized cells, they are very sensitive to their environment. Varied growth and differentiation properties can result in response to culture conditions (such as source and lot of FBS), time in culture (age related to passage number), confluency at time of passage and other environmental and laboratory factors, which can sub-select for cells with different growth characteristics<sup>3</sup>. For these reasons, optimizing the seeding density for the Caco-2 cells on the MultiScreen Caco-2 plate is highly recommended.

Laboratories that routinely use 24-well filter plates may elect to optimize seeding density for the 96-well filter plate based on the seeding density used for the 24-well plate. This can be accomplished by calculating the number of cells/cm<sup>2</sup> (see *Tips & Techniques*, at left, for converting cells/mL to cells/cm<sup>2</sup> to cells/well). Select routinely used seeding densities in cells/cm<sup>2</sup>, convert to cells/96-well filter area, and follow the methods described in this Protocol Note to evaluate the cells.

- a. Detach Caco-2 cells from cell culture flask using trypsin/EDTA as described in step **d** of *Section I, Medium, Additives and Caco-2 Cell Cultivation*, above.
- b. Once the cells are detached and resuspended in Caco-2 cell culture medium, count the cells to determine the cell number/mL.
- c. For a 10-day Caco-2 cell culture, divide cell suspension into 6 sterile, 15 mL centrifuge tubes. Dilute the cell aliquots with Caco-2 cell culture medium into cell concentrations ranging from 125,000 cells/mL to 500,000 cells/mL. These cell concentrations can be used to seed 9,375 cells/well to 37,500 cells/well or 85,000 cells/cm<sup>2</sup> to 341,000 cells/cm<sup>2</sup>. This is based on a seeding volume of 75 µL per filter well (see step **d** below). See *Tips & Techniques*, at left, for calculations and Figure 4, on page 5, for an example.


*Note: Achieving a uniform cell suspension when initially plating the cells will promote a more consistent monolayer across the 96 wells. This may be particularly difficult when seeding multiple plates. Frequent mixing is recommended to minimize the risk of large clumps of cells settling to the bottom of the tube, which could result in an inaccurate distribution of cells across the wells or plates.*

- d. Dispense 75 µL of cell dilution into the filter wells of the MultiScreen Caco-2 plates (96-well filter plate from Millipore cat. MACA C02 S5, MACA C02 S2 or MACA C02 B5). To optimize the seeding protocol, it is recommended to set up several columns for each cell density (minimally, 16 wells per seeding concentration). See Figure 4, on page 5, for an example.

- e. Dispense 250  $\mu\text{L}$  of Caco-2 cell culture medium into each of the 96 wells of the receiver plate (MACA C02 S5), or alternatively place 35 mL in the single-well feeder plate (MACA C02 B5). This may be accomplished by dispensing medium through the basolateral access holes (see Figure 3) for the receiver plate or the large access holes located at **A1** and **H12** for the single-well, feeder plate. Alternatively, disassemble the filter plate from the receiver plate or single-well, feeder plate. Place the filter plate on a sterile surface in a laminar flow hood, and add medium directly to filter plate and receiver plate. Gently reassemble the two components, and place in the cell culture incubator.
- f. For a 21-day, Caco-2 cell culture, divide cell suspension into 4 to 6 sterile 15 mL centrifuge tubes. Dilute the cell aliquots with Caco-2 cell culture medium into different cell concentrations ranging from 80,000 cells/mL to 167,000 cells/mL. (These cell concentrations can be used to seed 6,000 cells/well to 12,500 cells/well or 55,000 cells/cm<sup>2</sup> to 114,000 cells/cm<sup>2</sup>. This is based on a seeding volume of 75  $\mu\text{L}$  per filter well. See step **g** below). See *Tips & Techniques*, page 3, for calculations and Figure 5, page 5, for an example.
- g. Dispense 75  $\mu\text{L}$  into the filter well of the 96-well MultiScreen Caco-2 plates. For better optimization, set up several columns for each cell density (minimally 16 wells per seeding concentration).
- h. Dispense 250  $\mu\text{L}$  of Caco-2 cell culture medium into each of the 96 wells of the receiver plate, or alternatively, place 35 mL in the single-well, feeder plate as described in step **e**.
- i. Incubate plates for 10 or 21 days. Replace the medium in the filter well and the growth plate every other day, beginning no sooner than 48 hours after initial plating. Refer to *Section III, Exchanging the Medium in the MultiScreen Caco-2 Assay System*.

*Note: Cells seeded on the MultiScreen Caco-2 plate should be placed in an incubator that provides adequate humidity control. If a significant difference is observed in the performance of the wells on the perimeter of the plate relative to the interior, this indicates that the culturing environment is not adequately humidified. A cell culture incubator with electronic humidity control is recommended. If this is not possible, place the plates in an incubator that is opened infrequently.*

**Figure 3:** Patented Design of the MultiScreen Caco-2 Plate




### III. Exchanging the Medium in the MultiScreen Caco-2 Assay System


Long-term cell culture requires exchange of medium to (1) remove waste products that may accumulate and become toxic, and (2) to replenish nutrients. Cells grown on the MultiScreen Caco-2 plates for longer than 3 or 4 days typically need growth medium exchanged to sustain long-term viability of the cells. The MultiScreen Caco-2 assay system has been successfully evaluated for the growth and differentiation of Caco-2 cells from as little as 7 days to as long as 24 days. Medium should be exchanged every 48 to 72 hours beginning no earlier than 48 hours after initial plating. Exchange the medium in the filter well (apical compartment) and receiver plate or single-well feeder plate (basolateral compartment). Alternatively, medium can be exchanged in the basolateral compartment and replenished in the apical compartment (bring the volume back to approximately 75  $\mu$ L to replace evaporated medium lost during incubation) without removing the medium remaining in the filter well.

- a. The most critical part of removing and replacing the medium in the MultiScreen Caco-2 assay system is to avoid damage to the cell monolayer and the filter on which it is supported. In these experiments, a multichannel manifold (Sigma cat. M2656) was used to aspirate the medium, and an electronic multiwell pipettor with speed control at low or medium speed (Matrix Impact® pipettor, cat. 6002) was used to dispense the fresh culture medium.

**Figure 4: 10-day Caco-2 Cell Seeding Density Template**



**Figure 5: 21-day Caco-2 Cell Seeding Density Template**



## Tips & Techniques

### Caco-2 Cell Differentiation

1. If no significant expression and localization of P-glycoprotein (P) to the apical plasma membrane is observed after 10 days in culture of Caco-2 cells, be certain that the seeding density chosen is appropriate. Typically, to achieve good differentiation of the Caco-2 cells in 10 days, it is important to start with a sub-confluent flask of cells and to seed the wells at a higher density than would normally be used for a 21-day culture of Caco-2 cells.
2. The growth and differentiation of the Caco-2 cell lines need to be carefully monitored when optimizing the assay for use in a drug transport analysis. Many factors may contribute to the variability in the cell behavior. Cell passage number, culture medium and seeding density can all influence how the cells perform on the MultiScreen Caco-2 plates. Both tight junction formation and polarized expression of membrane proteins can exhibit changes as the cell passage number increases. How this will ultimately affect the rate measured in the drug transport analysis needs to be understood in each laboratory. As a rule, do not passage the Caco-2 cells more than 25 or 30 times before thawing a new vial.<sup>5</sup>

- b. The first option for feeding the plates consists of removing the filter plate (the filter plate can be placed directly on the cell culture hood surface) from the feeding plate (receiver plate or single-well plate) and aspirating the medium from both plates using the multichannel manifold. Be careful to guide the tips of the manifold down the lower left side of the filter wells along the apical assist (see Figure 3, on page 4). Replace the medium by gently pipetting 75  $\mu$ L into the filter well along the apical assist, using an electronic multichannel pipettor. Replace the medium in the basolateral feeding compartments, either 250  $\mu$ L/well for the 96-well receiver plate or 35 mL in the single-well feeder plate. Gently reassemble the filter plate with the receiver plate or the single-well feeder plate.

*Note: For left handed users, it may be more comfortable to rotate the plate 180° so that the apical assist is in the upper right corner.*

- c. The second option is to exchange the medium without disassembling the plate components. Aspirate the volume from the feeder plate directly through the basolateral access holes (see Figure 3, on page 4) or, if using the single-well feeder plate, aspirate through the large access holes located adjacent to either A1 or H12. Guide the tips of the aspirator into the holes (preferably using the recommended multichannel manifold), and carefully aspirate the medium from the filter wells. Take care not to contact the cell monolayer during aspiration.
- d. Replace the medium in the filter wells as directed in step b, and replace the medium in the feeding plate using the basolateral access holes for the receiver plate or the large access holes at either A1 or H12 for the single-well feeding plate.

## IV. Evaluating Caco-2 Cell Seeding Density

Once the Caco-2 cells have been in culture for the desired length of time, verify the integrity of the cell monolayer by measuring the transepithelial electrical resistance (TEER) for every well. This is a non-invasive method for determining monolayer integrity and can be used prior to the addition of test samples. Wells which have the appropriate electrical resistance can then be used as test wells for the transport studies. Another method for verifying monolayer integrity utilizes the fluorescent dye, lucifer yellow (LY) (Sigma cat. L0144). After the completion of the drug transport experiment, LY can be added to each well for %LY passage studies. Alternatively, a control population of cells within a plate can be selected for addition of the dye to monitor monolayer integrity during drug transport experiments.<sup>4</sup>

### Transepithelial Electrical Resistance Measurement

- a. At the end of the desired growth period, remove the plates from the incubator and allow them to equilibrate to room temperature, approximately 1 hour. Measure the electrical resistance across the monolayer using the Millicell®-ERS system ohm meter (Millipore cat. MERS 000 01) and the STX-100M electrode (World Precision Instruments, Sarasota, FL).
- b. Position the probe such that one end is immersed in the medium inside the filter well and the other is placed through the basolateral access hole into the medium in the growth plate. The probe should sit flat on the plate when positioned correctly with the thin collared electrode in the basolateral access hole. Record the electrical resistance for each well.

### Lucifer Yellow Passage Test

- a. Rinse the monolayer three times with 100  $\mu$ L HBSS (Invitrogen cat. 14025) using the same method described for *Section III, Exchanging the Medium in the MultiScreen Caco-2 Assay System*, substituting 100  $\mu$ L HBSS for the 75  $\mu$ L cell culture medium.
- b. Add 100  $\mu$ L of lucifer yellow at a concentration of 100  $\mu$ g/mL to each well in the filter plate after transport studies or to selected wells in a plate prior to transport studies.
- c. Add 250  $\mu$ L HBSS to the basolateral compartments of a 96-well transport analysis plate.
- d. Assemble the plate components and incubate for 1 hour at 37 °C.
- e. Remove the filter plate from the transport analysis plate and place into a fluorescent plate reader. [Data obtained in this Protocol Note were collected using a Wallac Victor<sup>2</sup>™ plate reader (PerkinElmer, Boston, MA).]
- f. Determine the fluorescent absorbance using wavelengths of 485 nm excitation and 535 nm emission.
- g. Calculate the percent of lucifer yellow passage across the cell monolayer into the transport analysis plate. Use the relative fluorescence measured for an equilibrium dilution of lucifer yellow in a separate analysis plate for reference. Following is the procedure:

#### Calculation of Percent of Lucifer Yellow Passage

Measure the relative fluorescence units (RFU) in 250  $\mu$ L of a 23  $\mu$ g/mL solution of lucifer yellow. This is the equilibrium dilution of the starting material. If the RFU measured in the basolateral compartment is equal to the RFU measured in this sample, 100% passage of lucifer yellow occurred:

$$\frac{(100 \mu\text{g/mL})(0.075 \text{ mL})}{0.075 \text{ mL} + 0.250 \text{ mL}} = 23 \mu\text{g/mL} = \text{RFU (equilibrium)}$$

In addition, determine the RFU for HBSS alone (blank). Then use these values to calculate the lucifer yellow passage in the test wells using the following equation:

$$\% \text{ Lucifer Yellow Passage} = \left[ \frac{\text{RFU (test)} - \text{RFU (blank)}}{\text{RFU (equilibrium)} - \text{RFU (blank)}} \right] \times 100$$

For example, if the measured values for each of these solutions equaled RFU (test) = 2000, RFU (blank) = 1000, and RFU (equilibrium) = 300,000, then the percent lucifer yellow passage would equal:

$$\left[ \frac{2000 - 1000}{300,000 - 1000} \right] \times 100 = 0.3\%$$

Choose the seeding density that provides the highest average electrical resistance with the least variability (lowest CV) and lowest lucifer yellow passage.



---

## V. Staining Cell Monolayers with Hema 3® Stain Kit

To visualize cell growth on the filter at the completion of the TEER and lucifer yellow testing, it may be useful to stain the cell monolayers to confirm the integrity and uniformity of the cells. One method for doing this is to use the Hema 3 stain kit (cat. 22122911) available from Fisher Scientific.

- a. Aspirate the wells to be stained. Add 100  $\mu$ L of fixative and incubate 1 minute at room temperature.
- b. Aspirate the fixative and add 100  $\mu$ L Solution 1. Incubate for 5 minutes at room temperature.
- c. Aspirate Solution 1 and add 100  $\mu$ L Solution 2. Incubate for 5 minutes at room temperature.
- d. Wash the wells extensively with Milli-Q® water, aspirate any residual liquid and allow wells to dry overnight.
- e. Visualize the wells using a stereomicroscope.

## VI. Performing a Drug Transport Assay using the MultiScreen Caco-2 Assay System

The ultimate goal for using the MultiScreen Caco-2 assay system is to perform a drug transport assay. This section provides a guide for the volumes to be tested in the device. In addition, washing the monolayer prior to the addition of test compounds is recommended. For drug transport experiments using radio-labeled drugs, unlabeled known drug (10  $\mu$ M) was spiked with a trace amount of  $^3$ H radio-labeled drug to monitor the drug transport rate. [Unlabeled mannitol (cat. M9546), digoxin (cat. D6003), propranolol (cat. P0884), testosterone (cat. T1500) and methotrexate (cat. M8407) were purchased from Sigma (St. Louis, MO). Radio-labeled mannitol (cat. NET101), digoxin (cat. NET222), propranolol (cat. NET515) and testosterone (cat. NET553) were purchased from PerkinElmer (Boston, MA). Radio-labeled methotrexate (cat. MT-701) was purchased from Moravek (Brea, CA).]

- a. When the 10-day or 21-day Caco-2 cultures have reached confluence and are differentiated, they are ready to be used for transport studies. Remove the MultiScreen Caco-2 plate from the incubator, and determine the electrical resistance for each well (as described in *Section IV, Evaluating Caco-2 Cell Seeding Density*). Next, wash the monolayer, exchanging the volume three times using sterile HBSS, pH 7.4. The method for the addition and removal of the wash buffer is similar to exchanging medium (*Section III, Exchanging the Medium in the MultiScreen Caco-2 Assay System*) except that the volume in the filter well may be increased to 100  $\mu$ L/well.
- b. Transfer the filter plate to a 96-well transport analysis plate after washing is complete.
- c. To determine the rate of drug transport in the apical to basolateral direction, add 75  $\mu$ L of the test compounds to the filter well. Drug concentrations typically ranging from 10  $\mu$ M to 200  $\mu$ M and diluted in HBSS, pH 7.4 or in an alternative buffer of desired pH, may be used. Fill the wells in the transport analysis plate with 250  $\mu$ L buffer.



- d. To determine transport rates in the basolateral to apical direction, add 250  $\mu\text{L}$  of the test compounds to the transport analysis plate wells. Fill the filter wells (apical compartment) with 75  $\mu\text{L}$  of the buffer.
- e. Incubate at 37 °C with or without shaking at 60 rpm on a rotary shaker. Typical incubation times are 1 to 2 hours. The data in this protocol note were generated using shaking for 1 hour.
- f. At the end of the transport period, disassemble the plates or remove a fixed volume (typically 50  $\mu\text{L}$ ) directly from the apical and basolateral wells (using the basolateral access holes) and transfer to a fresh transport analysis plate for LC/MS analysis.
- g. For radio-labeled drug evaluation, remove 20  $\mu\text{L}$  from each compartment and transfer to a plate (Wallac/Perkin Elmer 96-well flexible plate, cat. 1450-401) containing 100  $\mu\text{L}$  scintillation fluid. Mix and determine the radioactivity per sample using a multiwell plate scintillation reader such as the Wallac 1450 Microbeta® Plus Scintillation Counter or the Microbeta Trilux Multiwell Plate Scintillation Counter from PerkinElmer (Boston, MA).

## VII. Calculating Drug Transport Rates

The apparent permeability ( $P_{app}$ ), in units of centimeter per second, can be calculated for Caco-2 drug transport assays using the following equation:<sup>1</sup>

$$P_{app} = \left( \frac{V_A}{\text{Area} \times \text{time}} \right) \times \left( \frac{[\text{drug}]_{\text{acceptor}}}{[\text{drug}]_{\text{initial, donor}}} \right)$$

Where  $V_A$  is the volume (in mL) in the acceptor well, **Area** is the surface area of the membrane (0.11 cm<sup>2</sup> for MultiScreen Caco-2 plate), and **time** is the total transport time in seconds. For radio-labeled drug transport experiments the CPM units obtained from the Microbeta Trilux Multiwell Plate Scintillation Counter are used directly for the drug acceptor and initial concentrations such that the formula becomes:

$$P_{app} = \left( \frac{V_A}{\text{Area} \times \text{time}} \right) \times \left( \frac{\text{CPM}_{\text{acceptor}}}{\text{CPM}_{\text{initial, donor}}} \right)$$

## Results and Discussion

### Optimal Seeding Density for Caco-2 Cells

The optimal seeding density was determined for a 10- and 21-day culture of Caco-2 cells using the method described in this Protocol Note. Six different seeding densities were tested, with 32 wells tested per seeding density (see Figure 4, on page 4, and Figure 5, on page 5). The electrical resistance and lucifer yellow passage were determined for the 10- and 21-day cultures of Caco-2 cells (see Figure 6, below).

The highest mean electrical resistance (1700 and 1801 ohms) with the lowest variability (4.4 and 5.9 %CV) for the 10-day culture was observed at 15,000 and 20,000 cells/well. The lucifer yellow passage was lowest (0.3%) with the lowest variability at these same dilutions.


Therefore, seeding MultiScreen Caco-2 plates with Caco-2 cells for a 10-day culture would likely be successful starting with either 15,000 or 20,000 cells/well (which equals 200,000 to 267,000 cells/mL or 136,000 to 180,000 cells/cm<sup>2</sup>).

The highest mean electrical resistance (1999 and 2034 ohms) with the lowest variability (10 and 7.4 %CV) for the 21-day culture was observed at 9,000 and 11,000 cells/well. The lucifer yellow passage was lowest (0.2%) with the lowest variability at these same dilutions. Therefore, seeding MultiScreen Caco-2 plates with Caco-2 cells for a 21-day culture would likely be successful starting with 9,000 to 11,000 cells/well (which equals 120,000 to 147,000 cells/mL or 82,000 to 100,000 cells/cm<sup>2</sup>).

### Drug Transport Study Results

Caco-2 cells were cultivated for 10- or 21-days in the MultiScreen Caco-2 assay system using optimal seeding densities (20,000 and 9,000 cells/well, respectively), and a drug transport assay was performed as in Section VI. Drugs (mannitol, digoxin, propranolol, testosterone and methotrexate) at 10  $\mu$ M were added either apically, in HBSS supplemented with 10 mM MES, pH 6.8 or basolaterally, in HBSS at pH 7.4 and incubated for 1 hour at 37 °C, shaking at 60 rpm. The TEER was measured prior to the start of the drug transport assay, and lucifer yellow passage was evaluated on selected wells within the plate during the drug transport assay. The results demonstrate that the monolayers were well formed, the TEER values were higher at 21 than 10 days (as

**Figure 6:** TEER and Lucifer Yellow Passage Results for 10- and 21-day Caco-2 Optimization Experiments



**Table 1:** Performance of 10- and 21-day Caco-2 Cell Cultures at Optimized Seeding Densities on MultiScreen Caco-2 Plates

	TEER ohms-cm <sup>2</sup>	% Lucifer Yellow Passage	Drug Transport P <sub>app</sub> (10 <sup>-6</sup> cm/s) *				
			Mannitol	Propranolol	Digoxin	Testosterone	Methotrexate
10 day	158	0.5	1.05 ± 0.11	9.2 ± 1.0	1.1 ± 0.1	34.5 ± 2.3	0.6 ± 0.08
21 day	239	0.1	1.61 ± 0.24	8.9 ± 1.8	0.36 ± 0.09	33.4 ± 3.3	0.9 ± 0.16
Permeability			Low <sup>1</sup>	High <sup>1</sup>	Low <sup>6</sup>	High <sup>1</sup>	Low <sup>7</sup>

\*Note: The data shown are selected from three plates tested on the same day with a representative average from one plate shown. The standard deviation is calculated from the 12 replicates per drug.

expected), the lucifer yellow passage was less than 1%, and the drugs tested gave transport rates which provided appropriate classification into high and low permeability (see Table 1). In addition, digoxin is appropriately categorized as an active efflux drug since the digoxin transport in the basolateral to apical direction was shown to be 22 times faster than in the apical to basolateral direction after 10 days in culture. By 21 days, this difference increased to 60 times faster (data not shown).

The FDA Biopharmaceutical Classification system describes the use of the Caco-2 cell model for drug absorption measurements. The system characterizes new drug compounds as either low or high permeability. The drug transport rates measured here demonstrate that the MultiScreen Caco-2 assay system can be used to determine whether a drug is highly permeable by assigning low permeability to compounds with apparent permeability rates less than  $2 \times 10^{-6}$  cm/s.

The MultiScreen Caco-2 assay system is an excellent high throughput tool to screen compounds for intestinal transport properties. It provides a robust and reproducible format for conducting automated, high-throughput studies. Once optimized, growing cells in the MultiScreen Caco-2, 96-well, assay system provides an opportunity to measure drug transport rates of known and unknown compounds in a controlled and consistent manner. Other polarized epithelial cell model systems can be cultured in this device and used for discovery screening and ADME applications (see Millipore Protocol Note PC1061EN00: *Optimization of MDCK cell growth and differentiation for drug transport assay studies using a 96-well assay system*).

## References

1. Artursson, P., Karlsson, J., *Correlation between oral drug absorption in humans and apparent drug permeability coefficients in human intestinal epithelial (Caco-2) cells*. *Biochem. Biophys. Res. Comm.*, 1991. **175**:880–885.
2. Artursson, P., *Epithelial transport of drugs in cell culture. I: A model for studying the passive diffusion of drugs over intestinal absorptive (Caco-2) cells*. *J. Pharm. Sci.*, 1990. **79**:476–482.
3. Artursson, P., Palm, K., and Luthman, K., *Caco-2 monolayers in experimental and theoretical predictions of drug transport*. *Adv. Drug Deliv. Rev.*, 2001. **46**:27–43.
4. Bailey, C.A., Bryla, P., Malick, A.W., *The use of intestinal epithelial cell culture model, Caco-2, in pharmaceutical development*. *Adv. Drug Deliv. Rev.*, 1996. **22**:85–103.
5. Arena, A.A., Kamath, L. and Phillips, J.E., *Development of a Caco-2 drug transport assay in 96-well format*. Poster, American Association of Pharmaceutical Sciences, Denver, CO, 2001.
6. Stephens, R.H., O'Neill, C.A., Warhurst, A., Carlson, G.L., Rowland, M. and Warhurst, G., *Kinetic profiling of P-glycoprotein-mediated drug efflux and human intestinal epithelia*. *JPET*, 2001. **296**:584 – 591.
7. Bergstrom, C.A., Stratford, M., Lazorova, L., Avdeef, A., Luthman, A., and Artursson, P., *Absorption classification of oral drugs based on molecular surface properties*. *J. Med. Chem.*, 2003. **46**:558 – 570.

---

## Related Application and Protocol Notes

- PC1061EN00: Optimization of MDCK cell growth and differentiation for drug transport assay studies using a 96-well assay system
- AN1727EN00: Drug transport assays in a 96-well system: reproducibility and correlation to human absorption
- 

### To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office:

In the U.S. and Canada, call toll-free **1-800-MILLIPORE (1-800-645-5476)**

In the U.S., Canada and Puerto Rico, fax orders to **1-800-MILLIFX (1-800-645-5439)**

Internet: [www.millipore.com](http://www.millipore.com)

Technical Service:  
[www.millipore.com/techservice](http://www.millipore.com/techservice)

**pure e-commerce**

Now you can buy  
Millipore products  
online @

[www.millipore.com/purecommerce](http://www.millipore.com/purecommerce)

Millipore, MultiScreen, Millicell and Milli-Q are registered trademarks of Millipore Corporation.  
Lab-Line is a registered trademark of Lab-Line Instruments, Inc.

Wallac Victor<sup>2</sup> is a trademark of PerkinElmer or an affiliated company.

Microbeta is a trademark of Wallac Oy.

ATCC is a registered trademark of American Type Culture Collection.

Fisher Scientific and Hema 3 are registered trademarks of Fisher Scientific LLC.

Impact is a registered trademark of Matrix Technologies.

Lit. no. PC1060EN00 Rev.- 8/03 03-213

© 2003 Millipore Corporation, Billerica, MA U.S.A.  
All rights reserved.

### Millipore Worldwide

---

#### AUSTRALIA

Tel. 1 800 222 111  
or (02) 9888 8999  
Fax (02) 9878 0788

#### AUSTRIA

Tel. 0820 87 44 64  
Fax 0820 87 44 44

#### BALTIC COUNTRIES

Tel. +358 203 05 645  
Fax +358 203 05 644

#### BELGIUM

Tel. 070 225 645  
Fax 070 225 644

#### BRAZIL

Tel. (011) 5548-7011  
Fax (011) 5548-7923

#### CANADA

Tel. (800) 645-5476  
Fax (800) 645-5439

#### CHINA, PEOPLE'S REPUBLIC OF

Beijing:  
Tel. (86-10) 8519 1250  
(86-10) 6518 1058  
Fax (86-10) 8519 1255  
Guangzhou:  
Tel. (86-20) 8752 0187  
(86-20) 8752 0173  
Fax (86-20) 8752 0172  
Hong Kong:  
Tel. (852) 2803 9111  
Fax (852) 2513 0313  
Shanghai:  
Tel. (86-21) 5306 9100  
Fax (86-21) 5306 0838

#### CZECH REPUBLIC

Tel. 02 2051 3841  
Fax 02 2051 4298

#### DENMARK

Tel. 7010 5645  
Fax 7010 5644

#### EASTERN EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF

Tel. +33 3 88 38 9536  
Fax +33 3 88 38 9539

#### FINLAND

Tel. 0203 05 645  
Fax 0203 05 644

#### FRANCE

Tel. 0825 045 645  
Fax 0825 045 644

#### GERMANY

Tel. 01 805 045 645  
Fax 01 805 045 644

#### HUNGARY

Tel. 01 209 04 33  
Fax 01 209 02 95

#### INDIA

Tel. 080 839 46 57  
Fax 080 839 63 45

#### IRELAND

Tel. +44 1923 816375  
Fax +44 1923 818297

#### ITALY

Tel. 848 845 645  
Fax 848 845 644

#### JAPAN

Tel. (03) 5442-9711  
Fax (03) 5442-9736

#### KOREA

Tel. (822) 551-0230  
Fax (822) 551-0228

#### LUXEMBOURG

Tel. +32 2726 88 40  
Fax +32 2726 98 84

#### MALAYSIA

Tel. 03-7957-1322  
Fax 03-7957-1711

#### MEXICO

Tel. (55) 5576 9688  
Fax (55) 5576 8706

#### THE NETHERLANDS

Tel. 0900 7 645 645  
Fax 0900 7 645 644

#### NORWAY

Tel. 810 62 645  
Fax 810 62 644

#### POLAND

Tel. 22 669 12 25  
Fax 22 663 70 33

#### PORTUGAL

Tel. +34 91 728 39 60  
Fax +34 91 729 29 09

#### PUERTO RICO

Tel. (787) 273-8495  
Fax (787) 747-6553

#### SINGAPORE

Tel. 6842-1822  
Fax 6842-4988

#### SPAIN

Tel. 901 516 645  
Fax 902 011 644

#### SWEDEN

Tel. 0771 200 645  
Fax 0771 200 644

#### SWITZERLAND

Tel. 0848 645 645  
Fax 0848 645 644

#### TAIWAN

Tel. 886-2-2792-9333  
Fax 886-2-2792-6555

#### U.K.

Tel. 0870 900 46 45  
Fax 0870 900 46 44

#### U.S.A.

Tel. (800) 645-5476  
Fax (800) 645-5439

#### IN ALL OTHER COUNTRIES

Millipore International  
Tel. +1 781 533 8622  
Fax +1 781 533 8630

# MILLIPORE