# REPRODUCIBILITY BEGINS HERE

## Milli-Q® IX

7003/05/10/15

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### Production flow schematic

### Distribution flow schematic


- |  |  |  |  |
|--|--|--|--|
|  Tap feed water             |  IPAK Gard cartridge        |  RO membrane                  |  Elix EDI module                |
|  ech <sub>2</sub> o UV lamp |  Storage tank (25/50/100 L) |  Vent filter                  |  ech <sub>2</sub> o ASM UV lamp |
|  Strainer filter            |  E-POD dispenser            |  Millipak/Biopak final filter |  3-way valve                    |
|  Conductivity cell          |  Solenoid valve             |  Pressure sensor              |  Pump                           |
|  Thermistor                 |  Flow controller            |  Resistivity cell           |  |

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
# Purification unit and storage tank specifications

## Purification unit

## Storage tank

## Tubing and port requirements

| Item  | Description           |
|---|-----------------------|
| Feed water connection                               | ½" male NPT/BSP/GAZ   |
| Distance from feed water port                       | Maximum 5 m           |
| Distance from purification unit to E-POD® dispenser | Maximum 5 m           |
| Distance from purification unit to tank             | Maximum 5 m           |
| Power entry   | Connection IEC 13     |
| ON/OFF switch                                       | Available on the unit |
| Water sensor port                                   | Maximum 3.3 VDC       |
| Tank level adapter port                             | Maximum 5 VDC         |
| Ethernet port                                       | IEEE P802.3           |


## Electrical connections and specifications

| Item                    | Description               |
|-------------------------|---------------------------|
| Power source voltage    | 100~240 VAC ± 10%         |
| Power frequency         | 50/60 Hz ± 2 Hz           |
| Power used              | Maximum 350 VA            |
| Power cord length       | 2.5 m, plug: IEC13 female |
| Operational temperature | 4~40°C                    |
| Altitude                | Up to 3000 m              |
| Tank material           | High purity polyethylene  |

## Weights

| Milli-Q® IX system type |           | Dry weight | Shipping weight | Operating weight |
|-------------------------|-----------|------------|-----------------|------------------|
| Purification unit       | 7003/7005 | 19.5 kg    | 22.5 kg         | 23.4 kg          |
|                         | 7010      | 22.2 kg    | 25.4 kg         | 27.1 kg          |
|                         | 7015      | 22.5 kg    | 25.7 kg         | 27.4 kg          |
| Storage Tank            | 25 L      | 6.7 kg     | 8.5 kg          | 31.7 kg          |
|                         | 50 L      | 7.6 kg     | 10.6 kg         | 57.6 kg          |
|                         | 100 L     | 10.9 kg    | 12.8 kg         | 110.9 kg         |

# E-POD® dispenser specifications

## Tube and port requirements

| Parameter   | Description                              |
|---|--|
| Dispenser tubing length   | 0.92 m                                   |
| Distance from purification unit to POD dispenser                                  | Maximum 5 m                              |
| Distance between two POD dispensers<br>(Maximum 2 dispensers connected in series) | Maximum 5 m                              |
| POD dispenser data connection with unit   | Ethernet                                 |
| Electrical connection   | Powered by purification unit (24–28 VDC) |
| Foot pedal connection   | 3.3 V (5 m cable)                        |

## Weights

|        | Dry weight | Shipping weight | Operating weight |
|--------|------------|-----------------|------------------|
| E-POD® | 4.7 kg     | 7.2 kg          | 5.5 kg           |

## Screen description and functionalities

| Parameter              | Description   |
|------------------------|---|
| Capacitive touchscreen | Size: 5 in; Resolution: 800 × 480   |
| USB port               | USB 2.0 Highspeed standard  |
| Speaker                | Impedance: 8 Ω/max output power: 0.5 W  |
| Display in 9 languages | Chinese/English/French/German/Italian/Japanese/<br>Portuguese/Russian/Spanish |

# Water specifications

## Feed water requirements

|                                  |                                       |
|----------------------------------|---------------------------------------|
| Feed water                       | Potable tap water                     |
| Pressure                         | 1~6 bar                               |
| Temperature                      | 5–35°C                                |
| Conductivity                     | <2000 $\mu\text{S}/\text{cm}$ at 25°C |
| Dissolved CO <sub>2</sub>        | <30 ppm                               |
| Free chlorine                    | <3 ppm                                |
| Fouling Index                    | <10                                   |
| pH                               | 4–10                                  |
| Total Organic Carbon (TOC)       | <2 ppm                                |
| Langelier Saturation Index (LSI) | <0.3                                  |
| Hardness (as CaCO <sub>3</sub> ) | <300 ppm                              |
| Silica                           | <30 ppm                               |

## Pure, Type 2 water specifications<sup>1</sup>

|                                  |  |
|----------------------------------|--|
| Resistivity at 25°C <sup>2</sup> | >5 M $\Omega$ ·cm; typically 10–15 M $\Omega$ ·cm  |
| Conductivity at 25°C             | 0.2 $\mu\text{S}/\text{cm}$ ; typically 0.1 $\mu\text{S}/\text{cm}$  |
| TOC                              | $\leq$ 30 ppb  |
| Production flow rate             | 3 L/h (Milli-Q® IX 7003)<br>5 L/h (Milli-Q® IX 7005)<br>10 L/h (Milli-Q® IX 7010)<br>15 L/h (Milli-Q® IX 7015) |

## From an E-POD® dispenser with final filter, the following water quality specifications are achieved:<sup>1</sup>

|                                    |  |
|------------------------------------|--|
| Particulates <sup>3</sup>          | No particles with size >0.22 $\mu\text{m}$ |
| Bacteria <sup>4</sup>              | < 0.01 cfu/mL (< 10 cfu/L)                 |
| Pyrogens (endotoxins) <sup>5</sup> | <0.001 EU/mL                               |
| RNases <sup>6</sup>                | <1 pg/mL                                   |
| DNases <sup>6</sup>                | <5 pg/mL                                   |
| Proteases <sup>6</sup>             | <0.15 $\mu\text{g}/\text{mL}$              |
| Flow rate                          | Up to 2 L/min                              |

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# International regulatory requirements

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# Ordering information

| System components  | Catalog number |
|--|----------------|
| Milli-Q® IX 7003 Water Purification System (3 L/h production flow rate)  | ZIX7003T0C     |
| Milli-Q® IX 7005 Water Purification System (5 L/h production flow rate)  | ZIX7005T0C     |
| Milli-Q® IX 7010 Water Purification System (10 L/h production flow rate) | ZIX7010T0C     |
| Milli-Q® IX 7015 Water Purification System (15 L/h production flow rate) | ZIX7015T0C     |
| E-POD® Remote Dispenser  | ZIQEP0D00      |
| Milli-Q® IX 7003 System (3 L/h) with E-POD® Remote Dispenser             | ZIX7003P0C     |
| Milli-Q® IX 7005 System (5 L/h) with E-POD® Remote Dispenser             | ZIX7005P0C     |
| Milli-Q® IX 7010 System (10 L/h) with E-POD® Remote Dispenser            | ZIX7010P0C     |
| Milli-Q® IX 7015 System (15 L/h) with E-POD® Remote Dispenser            | ZIX7015P0C     |
| Milli-Q® IQ Storage Tank Frame 25 L                                      | TANKA025       |
| Milli-Q® IQ Storage Tank Frame 50 L                                      | TANKA050       |
| Milli-Q® IQ Storage Tank Frame 100 L                                     | TANKA100       |
| Milli-Q® IQ Storage Tank Top Assembly (includes ASM)                     | TANKT0PA1      |

| Purification kits & Application POD-Paks                                     | Catalog number |
|--|----------------|
| Milli-Q® IX 7003/05 Pure Water Purification Kit (pretreatment & vent filter) | IX700XPKIT     |
| Milli-Q® IX 7010/15 Pure Water Purification Kit (pretreatment & vent filter) | IX70XXPKIT     |
| IPAK Gard® 03/05 Pretreatment Cartridge for Hard Water*                      | IPAKGARDH1     |
| IPAK Gard® 10/15 Pretreatment Cartridge for Hard Water*                      | IPAKGARDH2     |
| Vent Filter*   | TANKV01A1      |
| Vent Filter HF (for high-flow applications)* †                               | TANKVH1A1      |
| Millipak® 0.22 µm Filter   | MPGP002A1      |
| Millipak® Gold 0.22 µm Sterile Filter  | MPGPG02A1      |
| Biopak® Polisher   | CDUFBIOA1      |

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| 액세서리  | Catalog number |
|---|----------------|
| Foot pedal  | ZMQSFTSA1      |
| Sanitary sampling valve kit                       | ZIQ7ESP01      |
| System wall mounting bracket                      | SYSTFIXA1      |
| E-POD® wall mounting bracket                      | WMBQP0D01      |
| Tank wall mounting bracket                        | TANKFIXA1      |
| Water sensor                                      | ZWATSENA1      |
| External solenoid valve for feed water            | EXTSV00A1      |
| Washer distribution kit 230 V (right)             | ZWDK5R100      |
| Washer distribution kit 230 V (left)              | ZWDK5L100      |
| Washer distribution kit 115 V (right)             | ZWDK6R100      |
| Washer distribution kit 115 V (left)              | ZWDK6L100      |
| Wall mounting bracket for washer distribution kit | WMBWASH1       |

| System care                         | Catalog number |
|-------------------------------------|----------------|
| ROCare A - Acidic care              | ZWACID012      |
| ROCare B - Basic care               | ZWBASE012      |
| ROProtect C - Chlorine tablets      | ZWCL01F50      |
| EfferSan Effervescent Tablets (USA) | 5874316024     |
| EfferSan Effervescent Tablets (CAN) | 5874316024C    |

| Connections                   | Catalog number |
|-------------------------------|----------------|
| Connector 2m system-to-POD    | ZFC0NN2SQ      |
| Connector 5m system-to-POD    | ZFC0NN5SQ      |
| Connector 2m system-to-tank   | ZFC0NN2ST      |
| Connector 5m system-to-tank   | ZFC0NN5ST      |
| Connector 2m(POD-POD)         | ZFC0NN2QQ      |
| Connector 5m POD-to-POD       | ZFC0NN5QQ      |
| Multi-system installation kit | ZIQ7MSKT1      |

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