3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**

# **MS RT Calibration Mix**

Catalog Number MSRT1 Storage Temperature -20 °C

# **Product Description**

Standards are critical in mass spectrometry-based proteomics to ensure optimal and consistent system performance before, during, and after sample analysis. They may be used to assess peptide elution, troubleshoot chromatography, predict retention time, and to demonstrate that the LC-MS platforms are working properly.

MSRT1 is a mixture of 14 isotopically labeled synthetic peptides designed to act as a LC-MS platform standard, containing elements to test:

- LC Resolution •
- Peptide elution profile •
- LC-MS Platform performance •
- -Retention time prediction

The peptides in MSRT1 are designed to span the normal elution profile of complex proteomic samples. Also, the peptides have been added in various amounts (based on electrospray response) to have relatively similar responses and be readily observed.

MSRT1 can be used as a guick test of LC-MS platform performance, comparison of LC gradients and columns, monitoring column and system changes, or used to predict retention times of peptides across LC-MS platforms.

# Reagent

MS RT Calibration Mix

1 vial

Each vial contains various amounts of 14 isotopically labeled synthetic peptides dried via vacuum centrifugation. Each peptide is isotopically labeled with  $({}^{13}C, {}^{15}N)$  Leucine (+7 Da),  $({}^{13}C, {}^{15}N)$  Lysine (+8 Da), or  $({}^{13}C, {}^{15}N)$  Arginine (+10 Da).

# Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

# Storage/Stability

The dried peptide product is stable for up to 2 years when stored at -20 °C.

# **Preparation Instructions**

The peptides are soluble in reagents compatible with proteomic workflows (e.g., 0.1-1.0% formic acid). However peptide recovery is maximized when ~20% organic solvent is used to resuspend the dried material. Reconstitution with <20% CH<sub>3</sub>CN may reduce recovery of certain peptides.

# Procedure

- 1. Briefly centrifuge the vial at  $\sim 10,000 \times g$  to collect the product at the bottom.
- 2. Add  $\geq 10 \ \mu L$  of 20% CH<sub>3</sub>CN and 1.0% formic acid.
- Vortex rigorously for 30–60 seconds.
- Briefly centrifuge the vial at ~10,000  $\times$  g to collect 4. the product at the bottom.
- 5. Add at least  $10 \times$  the volume used in step 2 of 0.1% formic acid. <u>Note</u>: If the total volume is 100  $\mu$ L, then the concentration (based on the amount of GLFIIDD[K] peptide) is 1 pmol/ $\mu$ L. This is the recommended storage concentration. The solution can be stored up to 1 month at 2-8 °C. DO NOT FREEZE after reconstitution.
- 6. Vortex rigorously for 30-60 seconds.
- 7. Dilute desired amount of material prepared in step 5 in 0.1% formic acid as suitable for the particular LC-MS analytical system. Note: It is recommended that steps 6 and 7 be performed on the day of analysis. Some loss of the hydrophobic peptides has been observed upon storage below 1 pmol/µL (based on the amount of GLFIIDD[K] peptide).


Instrument parameters used for selected reaction monitoring analysis (SRM) are provided via a Skyline file and an Excel table, available for download at www.sigmaaldrich.com. Users may optimize their specific SRM parameters as necessary.


# Table 1. MSRT1 Peptides

Peptide #	Peptide Sequence	Amount (pmols)	Monoisotopic MW	iRT *
1	RGDSPASSP[K]	132	1008.5080	-29.6
2	GLV[K]	309	423.2937	-27.4
3	LGGNETQV[R]	81	982.5071	-25.0
4	AEFAEVS[K]	101	887.4480	-9.4
5	SGFSSVSVS[R]	67	1021.5068	4.3
6	ADEGISF[R]	76	903.4325	9.0
7	DISLSDY[K]	120	947.4691	21.0
8	LVNEVTEFA[K]	38	1156.6219	31.4
9	DQGGELLSL[R]	37	1096.5752	38.8
10	GLFIIDD[K]	100	927.5157	54.6
11	LGEYGFQNA[L]	82	1117.5517	61.4
12	YWGVASFLQ[K]	61	1205.6324	84.3
13	TDELFQIEGLKEELAYL[R]	350	2176.1291	111.4
14	AVQQPDGLAVLGIFL[K]	46	1675.9752	126.7

Amino acid in [brackets] denotes site of label incorporation as follows: [K],  ${}^{13}C_{6}{}^{15}N_{2}$ ; ; [R],  ${}^{13}C_{6}{}^{15}N_{4}$ ; or [L],  ${}^{13}C_{6}{}^{15}N_{1}$ . \* As determined by C18 chromatography followed by ESI mass spectrometry. iRT values generated using Skyline.<sup>1</sup>

#### Figure 1.





Dual Column Ascentis Express Peptide ES C18 (1 mm I.D.  $\times$  15 cm length, 2.7  $\mu$ m particles, 40 °C) at 90  $\mu$ L/min Solvent A: 0.1% formic acid Solvent B: Acetonitrile with 0.1% formic acid

#### Reference

1. MacLean, B., *et al.*, Skyline: An Open Source Document Editor for Creating and Analyzing Targeted Proteomics Experiments. *Bioinformatics*, **26**(7), 966-968 (2010).

This product is licensed under U.S. Patent No. 7,396,688 and foreign counterparts from E. I. du Pont de Nemours and Company. The purchase of this product conveys to the buyer the nontransferable right to use the purchased amount of the product for research and development only, including services for a third party for consideration. The buyer cannot sell or otherwise transfer this product, its components or materials made using this product or its components to a third party. Information about licenses for excluded uses is available from: E. I. du Pont de Nemours and Company; Attn: Associate Director, Commercial Development; DuPont Experimental Station E268; 200 Powdermill Rd.; Wilmington, DE 19803; 1-877-881-9787 (voice), 1-302-695-1437 (fax), licensing@dupont.com.

JJW ,PM,GCY,KR,IC,MAM 08/15-1

©2015 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.