

Preparing for the Paradigm Shift in Bioprocessing

While upstream titers have improved dramatically over the last decade, downstream processing has remained relatively unchanged. Now is the time to push bioprocess efficiencies even further. In August 2017, Merck acquired Natrix Separations, Inc., a provider of hydrogel membrane products for single-use chromatography. Renaud Jacquemart has been with Natrix since 2011 and today is the Director of Vaccines Process Sciences. Here, he explains how the Natrix® technology could offer a helping hand in next generation bioprocessing.

By Renaud Jacquemart

The exact definition of “next generation bioprocessing” and how that might be implemented varies depending on which company you speak to, but it essentially comes down to improving processes. And I believe it represents a paradigm shift in the industry. Biologic products are not simple to make and biopharma is a conservative industry, so there is not always a willingness to adopt new technologies. However, upstream titers have improved so much that many companies are now looking to embrace change; for example, by retrofitting their facilities to cope with the increased quantities of antibody that they are producing upstream. Some are producing kilograms of antibody per batch, whereas ten years ago they were only producing a few hundred grams; most downstream processes simply aren't built



to handle such large quantities, which leads to bottlenecks in purification. In addition to process challenges, there is also increasing pressure on drug costs, with governments, insurance companies and other payers calling for biopharma companies to lower their prices. One example that highlights the issue is Ocrelizumab (Ocrevus, Genentech), approved in May 2017 by the FDA to treat primary progressive multiple sclerosis, with a cost of \$65,000 per year. Such high prices put tremendous pressure on health systems, and therefore treatments are discussed not in terms of disease-modifying performance, but in improvement over cost ratios. My team and I have no direct impact on the pricing process but we can at least help manufacturers to develop the best processes possible, which may enable them to lower costs. After all, if you want to sell something at a lower cost, it stands to reason that you first have to produce it at a lower cost.

Continuous manufacturing is a buzzword in the biopharma industry right now because of its potential to intensify processes. By compressing the process into a small footprint that still produces the same quantities as a large-scale facilitate, you benefit from lower

upfront investment, lower running costs, and a shorter build time. Right now, however, there are still many questions about how to actually implement a fully continuous bioprocess. Single use is certainly a key enabler of continuous bioprocessing, and in upstream there has been a lot of progress with the introduction of perfusion bioreactors. There are also single-use options for certain unit operations in the downstream process, but one of the missing links was chromatography, which is one of the most expensive steps of the process. In the current paradigm, traditional resin chromatography columns are oversized to match upstream productivity. They are very expensive and require amortization over many production lots, as well as cleaning, storage, validation and significant labor to operate. And they also omit flexibility and agility in the manufacturing.

Forging the way with technology So how does the Natrix® technology fit into next generation bioprocessing? Merck is committed to establishing a portfolio of next generation bioprocessing technologies that enable faster and more efficient biomanufacturing – and one

of these technologies is the Natrix® Membrane. The Natrix® platform is not the only chromatography option on the market, but our product is unique in that it offers extremely high productivity. Indeed, where other chromatography technologies provide either high binding capacity or high speed, the Natrix® technology features a very high density of binding sites throughout a macroporous 3D structure that provide high productivity, simultaneously. Effectively, you can achieve production-scale capacity – and maintain quality – in a much smaller footprint, allowing you to move away from traditional columns. The technology, which originally came from McMaster University in Ontario (Canada), is made using a single-step polymerization that provides lot-to-lot reproducibility to guarantee consistent performance. And because our technology is single use and fully validated for GMP use, it is treated as a consumable, ready to use in a plug-and-play format.

Before I joined Natrix, I had actually worked with other membrane chromatographic products. Back in the early 2000, when fears of an avian flu pandemic were rife, membrane chromatography was considered by many to be the holy grail because of its potential to accelerate vaccine development, but it didn't take off at the time because we didn't have the right capacities and the industry was not ready to fully adopt single-use technologies. I was really disappointed at the time, but today things are very different – intensified upstream processes deliver large quantities of biologics and patients are waiting; it is time to address the downstream bottlenecks!

The acquisition by Merck has greatly increased the visibility and attractiveness of our technology. Natrix Separations, Inc. is not a small company anymore – it is part of a large, well-established organization that is a leader in purification. I recently gave a presentation at an international, scientific

And Coming Up Next...

The work at Natrix is by no means complete – we are currently developing the second generation of our membranes that will benefit from the addition of affinity ligands. The first product we are working on is a protein A membrane. Protein A ligands have been designed to selectively capture antibodies, therefore reducing the number of process steps required and providing robustness to the process. But many other ligands and applications have largely been ignored. In the future, we envision that all end users will have access to a Natrix® affinity membrane column and achieve the same performance that is obtained with Protein A today. For example, it is possible to make vaccines with the same approach. In this case, a ligand is attached to the base membrane to capture the virus of interest. My team recently published on this strategy (1,2). The performance of the Natrix® products in this context is a huge advantage, but major improvements would not be possible without a holistic approach, combining progress in upstream and downstream, but also facilitating design and analytical/release capabilities.

We are concretely playing a role in helping to make biotherapeutics as well as vaccines more affordable.

conference and people congratulated me on the deal with Merck, with one individual telling me that they had been investigating the Natrix® technology for some time, but were concerned about the risks of us being a smaller company. Now that we are part of a bigger company, they are happy to move forward because we have provided them with the real solution to address their challenges in therapeutic supplies!

Despite biopharma being a conservative industry, I am seeing a lot of positive

We currently partner with two world leaders in biomanufacturing innovation: Batavia Biosciences (Leiden, Netherlands) and Univercells (Gosselies, Belgium), with whom we were awarded the Grand Challenge “Innovations in Vaccine Manufacturing for World Markets” by the Bill & Melinda Gates Foundation. In this project, the aim is to develop a new vaccine manufacturing platform that intensifies the process and cuts the cost down to 15 cents per dose. We will use our chromatography membrane platform, Univercells' process intensification and integration technologies, and Batavia's vaccine development and manufacturing capabilities. The initial target is to establish a microfacility for inactivated polio vaccine, but the platform concept could be applied to any viral vaccine to enable the production and distribution of high quality vaccines in low income countries.

References

1. M Vandersluis et al., “Achieving Intensification and Flexibility in Virus Purification with Next-Generation Chromatography Tools,” *Pharmaceutical Technology*, 17 39-44 (2017).
2. M Zhao et al., “Affinity Chromatography For Vaccines Manufacturing: Finally Ready For Prime Time?” *Vaccines* (in press).

activity that implies a willingness to embrace next generation bioprocessing. I am editing a book right now called *New Paradigms in Biomanufacturing*, but the outline keeps changing because my partners worldwide are making progress so fast! It is really encouraging and shows how much momentum there is around this topic.

Renaud Jacquemart is Director Vaccines Process Sciences at Merck, Burlington, Ontario, Canada.