

sigma-aldrich.com

3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

Percoll[®]

Catalog Numbers P1644, P4937, and P7828

Percoll[®] PLUS Catalog Numbers E0414

Storage Temperature 2-8 °C

Product Description

Percoll[®] is a classic medium for density gradient centrifugation of cells, viruses, and subcellular particles. Percoll consists of colloidal silica particles of 15-30 nm diameter (23% w/w in water), which have been coated with polyvinylpyrrolidone (PVP). The PVP coating renders the product completely non-toxic and ideal for use with biological materials. The PVP is firmly bound to the silica particles as a monomolecular layer. Due to its heterogeneity in particle size, sedimentation occurs at different rates, spontaneously creating very smooth, isometric gradients in the range of 1.0–1.3 g/ml. Most biological particles having sedimentation coefficient values greater than 60S can be successfully isolated in Percoll gradients.

Percoll PLUS is also a silica-based colloidal medium for cell separation by density gradient centrifugation. The silica particles of this medium are covalently coated with silane, providing product stability and long shelf life. The silane coating also provides low osmolality and toxicity, as well as low viscosity. Percoll PLUS has low levels of endotoxins, making it well-suited for cell separation in clinical research applications. Cell separation is performed using the same conditions as the original Percoll. The low toxicity of Percoll and Percoll PLUS ensures that removal of the medium from separated cellular particles is not usually necessary.

Figure 1.

Properties of Percoll/Percoll PLUS

Percoll PLUS/Percoll have the following features:

- Low osmolality permitting precise adjustment to physiological conditions without significant interference from the medium.
- Compatibility with living cells and viruses, allowing separation and recovery of intact, fully active systems.
- Impermeable to biological membranes, resulting in no change of buoyant density of particles during centrifugation.
- Spontaneous formation of gradient during centrifugation, allowing mixing of large sample volumes in the centrifuge tubes.
- Low viscosity resulting in rapid formation of gradients and particle separation.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Property	Percoll PLUS	Percoll
Composition	Silica sol with covalently linked silane	Silica sol with non-dialyzable PVP coating
Density (g/ml)	1.130 ± 0.005	1.130 ± 0.005
Osmolality (mOsm/kg H ₂ O)	maximum 30	maximum 25
Conductivity (mS/m)	-	maximum 100
Viscosity (cP)	maximum 15 (20 °C)	maximum 15 (20 °C)
рН	9.4 ± 0.5 (20 °C)	9.0 ± 0.5 (20 °C)
Endotoxin (EU/ml)	<2	-
(Supplier information)		

(Supplier information)

Preparation Instructions

Percoll/Percoll PLUS is best used in balanced salt solutions, physiological saline, or 0.25 M sucrose. Cells can be separated in gradients in balanced salts solutions. Subcellular particles, however, tend to aggregate in the presence of salts and it is recommended the separation of such particles be carried out in Percoll/Percoll PLUS diluted with sucrose (0.25 M final concentration).

The low osmolality of Percoll/Percoll PLUS permits this parameter to be controlled by the user without interference from the density medium itself. The addition of 9 parts (v/v) of Percoll/Percoll PLUS to one part (v/v) of either 1.5 M NaCl, 10× concentrated culture medium, or 2.5 M sucrose will result in a solution adjusted to ~340 mOsm/kg H₂O. Final adjustments can be made with the addition of salts or distilled water. The precise osmolality should be checked prior to use with an osmometer.

Percoll/Percoll PLUS can be used within the pH range of 5.5–10.0 without any changes in properties. Percoll/Percoll PLUS may form a gel at pH values below 5.5. Gelling can also be caused by the presence of divalent cations, particularly at elevated temperatures.

Percoll/Percoll PLUS will form self-generated gradients by centrifugation at $10,000 \times g$ in 0.15 M saline or $25,000 \times g$ in 0.25 M sucrose in fixed angle rotors after 15 minutes. Cells or subcellular particles can be mixed with Percoll/Percoll PLUS prior to centrifugation and will band isopycnically as the gradient is formed *in situ*. The use of swinging bucket rotors is not recommended for self-generating gradients.

Percoll/Percoll PLUS may be diluted directly to make a final working solution of known density by the following procedure. In a graduated cylinder, add 1.5 M NaCl or 2.5 M sucrose to 1/10 the desired volume. To this add the required calculated volume of undiluted Percoll/Percoll PLUS (see Figure 2). Make up to the final volume with distilled water.

Figure 2.

Volume Calculation

$$V_0 = V \times \rho - 0.1 \rho_{10} - 0.9 \rho_0 - 1$$

- V_o = Volume of undiluted Percoll/Percoll PLUS required in ml
- V = Volume of final working solution in ml
- ρ = Desired density of final working solution
- ρ_o = Density of Percoll/Percoll PLUS undiluted (lot specific value)
- $\label{eq:rho10} \begin{array}{l} \rho_{10} \mbox{ = Density of 1.5 M NaCl (1.058 g/ml) or 2.5 M} \\ \mbox{ sucrose (1.316 g/ml)} \end{array}$

The formula is useful for achieving densities that will be very close to the actual densities required. However, slight variations in densities and volumes may affect final density. For highly accurate density requirements, it is recommended to check and adjust the final density using a densitometer or refractometer. The refractive index of diluted Percoll/Percoll PLUS solutions has a linear correlation with the solution density (see Figure 3).




Figure 3. Refractive index as a function of density of a Percoll gradient (work from GE Healthcare Bio-Sciences AB, Uppsala Sweden)

Storage/Stability

Percoll/Percoll PLUS is aseptically filled and can be stored for up to two years in an unopened container. If stored at -20 °C, gradients form upon thawing, necessitating mixing the contents of the bottle before use.

Preformed gradients can be stored for weeks without a change in gradient shape, provided the gradient is kept under aseptic conditions and not physically disturbed.

Percoll/Percoll PLUS may be autoclaved at 120 °C for 30 minutes without any change in properties. Autoclaving of solutions must be carried out without addition of salts or sucrose. The presence of salts will cause Percoll/Percoll PLUS to gel and the presence of sucrose will cause caramelization.

Minimal contact with air should be maintained during autoclaving to avoid formation of solid particles at the Percoll/Percoll PLUS/air interface. This can be accomplished by using a narrow-top bottle. If particles do form, they may be removed by filtration or low speed centrifugation. If any significant evaporation occurs during autoclaving, the volume should be replenished with sterile water so that the density is not affected. The plastic bottles in which Percoll/Percoll PLUS is packaged are not autoclavable.

Procedures

Examples of Separations in Percoll/Percoll PLUS

Source	Density (g/ml)	Centrifugation Conditions
Rat Liver Cells		
Hepatocytes	1.07–1.10	30,000 × <i>g</i> (30 min)
Kupffer cells	1.05–1.06	30,000 × <i>g</i> (30 min)
Human Cells		
Thrombocytes	1.04–1.06	*
Lymphocytes	1.06–1.08	*
Granulocytes	1.08–1.09	*
Erythrocytes	1.09–1.10	*
E. coli	1.13	30,000 × <i>g</i> (20 min)
Virus		
Tobacco mosaic	1.06	100,000× <i>g</i> (45 min)
Equine abortion	1.08	40,000 × <i>g</i> (45 min)
Influenza	1.06	25,000 × <i>g</i> (25 min)
Organelles		
Mitochondria	1.09–1.11	50,000 × <i>g</i> (45 min)
Lysozomes	1.04–1.07	50,000 × <i>g</i> (45 min)
	1.08–1.11	50,000 × <i>g</i> (45 min)
Peroxisomes	1.05–1.07	63,000 × <i>g</i> (30 min)
Synaptosomes	1.04–1.06	50,000 × g (45 min)
Nuclei	1.08–1.12	100,000 × <i>g</i> (60 min)

Separation of blood cells is best carried out by preforming the gradient (starting density 1.09 g/ml) by centrifugation at $20,000 \times g$ for 20 minutes, then layering blood on top of the gradient. Then centrifuge at $1,000 \times g$ for 5 minutes in a swinging-bucket rotor, leaving the thrombocytes in the serum layer above the gradient; the serum layer can be removed with a pipette (rate-zonal separation). A further spin for 20 minutes at

 $1,000 \times g$ separates the other cell types at their isopycnic densities.

After centrifugation, the gradient fractions can be collected by puncturing a hole in the bottom of the tube. Another simple and convenient method is to collect the fractions from the tube by displacement with a dense medium such as undiluted Percoll/Percoll PLUS or a 60–65% sucrose solution.

Percoll/Percoll PLUS does not interfere with fluorescent activated cell sorting or with electronic cell counting instruments.

<u>Removal of Percoll/Percoll PLUS from cells</u> Living cells can be separated from Percoll/Percoll PLUS medium by washing with physiological saline (5 volumes saline to 1 volume of cell suspension). The washing may be repeated two to three times and the cells collected between each washing step by centrifugation at $200 \times g$ for 2–10 minutes.

For viruses and subcellular particles, which are too small to be pelleted by low speed centrifugation, the particles can be separated from Percoll/Percoll PLUS by high-speed centrifugation. The undiluted fraction is centrifuged at $100,000 \times g$ for two hours in a swinging-bucket rotor or 90 minutes in an angle-head rotor. The biological material remains above the hard pellet of Percoll/Percoll PLUS.

References

Information on physical properties and applications was obtained from the supplier.

Lymphocyte Separation

Stibenz, D. and Buhrer, C., Scand. J. Immunol., 39, 59-63, 1994.

Pistoia, V. et al, Stem Cells, 11, 150-155 (1993).

Giddings, J.C. et al., *Clin. Lab. Hematol.*, 2, 121-128 (1980).

Monocytes

Przepiorka, D et al., *Am. J. Clin. Pathol.*, 95, 201-206 (1991).

Osipovich, O.A. et al., *In vitro. Biull. Eskp. Biol. Med,* 113, 638-640 (1992).

Giddings, J.C. et al., *Clin. Lab. Hematol.*, 2, 121-128 (1980).

Erythrocytes

Pascual, M. et al., *Eur. J. Immunol.*, 24, 702-708 (1994).

Vanden Berg, J.J.M. et al., *Arch. Biochem. Biophys.*, 298, 651-657 (1992).

Rennie, C. et al., Clin. Chim. Acta, 98, 119-125 (1979).

Natural Killer Cells Warren, H.S. and Skipsey, L.J., Immunol., 74, 78-85 (1991). Saksela, E. et al., Immunological Rev., 71-123 (1979). Krishnaraj, R. Cell. Immunol., 141, 306-322 (1992). Neutrophils Arnould, T. et al., *Blood*, 83, 3705-3716 (1994). Read, R.A. et al., Surgery, 114, 308-313 (1993). Conway, E.M. et al., Blood, 80,1254-1263 (1992). Eosinophils Schweizer, R.C. et al., Blood, 83, 3697-3704 (1994). Burgers, J.A. et al., Blood, 81, 49-55 (1993). Blom, M., et al., Blood, 83, 2978-2984 (1994). Kupffer Cells Page, D.T. and Garvey, J.S., J. Immunol. Methods, 27, 159-173 (1979). Hepatocytes Dou, M. et al., Cryobiol., 29, 454-469 (1992). Page, D.T. and Garvey, J.S., J. Immunol. Methods, 27, 159-173 (1979). Obrink, B. et al., Biochem. Biophys. Res. Comm., 77, 665-670 (1977). Basophils Kepley, C. et al., J. Immunol. Meth., 175, 1-9 (1994). Arock, M. et al., Int. Arch. Allergy Immunol., 102, 107-111 (1993). Tanimoto, Y. et al., Clin. Exper. Allergy, 22, 1020-1025 (1992). Levdig Cells Syed, V. et al., in vitro. J. Endocrinol. Invest., 14, 93-97 (1991).Schleicher, R.L. et al., Biol. Reprod., 48, 313-324 (1993). <u>Spermatozoa</u> Chen. Y. et al., Int. J. Fertil., 37, 315-319 (1992). Bongso, A. et al., Archiv. Androl., 31, 223-230 (1993). Morales, P. et al., Human Reprod., 6, 401-404 (1991). **Bone Marrow** Avraham, H. et al., *Blood*, 79, 365-371 (1992). Genot, E. et al., Blood, 80, 2060-2065 (1992). Louache, F. et al., Blood, 78, 1697-1705 (1991). Macrophages Calhoun, W.J., J. Lab. Clin. Med., 117, 443-452 (1991). Calhoun, W.J. et al., Am. Rev. Respir. Dis., 145, 317-325 (1992). Narahara, H., et al., Endocrinol. Metab, 77, 1258-1262 (1993).

Mast Cells Kulmburg, P.A. et al., Exp. Med., 176, 1773-1778 (1992).Inagaki, N. et al., Life Sci., 54, 1403-1409 (1994). Kurosawa, M. et al., Int. Arch. Allergy Immunol., 97, 226-228 (1992). Thymocytes Hoshino. J. et al., Biochem. International, 27, 105-116 (1992). Fearnhead, H.O. et al., Biochem. Pharmacol., 48, 1073-1079 (1994). Sun, X.M. et al., Biochem. Pharmacol., 44, 2131-2137 (1992). Cohen, G.M. et al., J Immunol., 153, 507-516 (1994). Pancreatic Islets Buitrago, A. et al., Biochem. Biophys. Res. Commun., 79, 823-828 (1977). Endothelial Sbarbati, R. et al., Blood, 77, 764-769 (1991). Neurons Soldenberg, S.S., and De Boni, U., J. Neurobiol., 14, 195-206 (1983). Platelet Membranes Perret, B. et al., Biochim. Biophys. Acta, 556, 434-446 (1979). Hepatocyte Membranes Obrink, B. et al., Biochem. Biophys. Res. Commun., 77, 665-670 (1977). CHO Membranes Cezanne, L. et al., Biochim. Biophys. Acta, 1112, 205-214 (1992). Lysozomes Lindley, E.R. and Pisoni, R.L., Biochem. J., 290, 457-462 (1993). Kominami, E. et al., J. Biochem., 111, 278-282 (1992). Mitochondria Chemnitus, J.M. et al., Int. J. Biochem., 4, 589-596 (1993). Lopez-Mediavilla, C. et al., Exp. Cell Res., 203, 134-140 (1992). Granules Kjeldson, L. et al., Blood, 83, 1640-1649 (1994). Sengelov, H. et al., Biochem. J., 299,473-479 (1994). Nuclei Hahn, C. and Covault., J., Anal. Biochem., 190, 193-197 (1990).

Percoll is a registered trademark of GE Healthcare.

VC,PHC 08/15-1

©2015 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.