



NovaSeptum® Pre-loaded Disposable Holder

Product Descriptions

The NovaSeptum® pre-loaded disposable holder is available with various standard configurations using NovaSeptum® high purity bags, bottles or accurate volume (AV) units. Other configurations are available upon request. "For an optimal aseptic solution, use the NovaSeptum® pre-loaded disposable holder together with a 1.5" or 2" NovAseptic® connector. A certificate of quality is included in each box.



NovaSeptum® pre-loaded 1.5" Single Use holder

Specifications

NovaSeptum® Single-Use holder, TC 1.5 / 2", 2 mm

Holder and Magazine	PEI Black
Locking ring	PEI Orange
Screw	A2-A4
Environmental	
Design Temperature, TC holder	Steam in place cycle at up to 134 °C (273 °F) or Autoclave during 1 hour at up to 134 °C (273 °F) (when applicable)
Traceability	The product and packaging label includes the catalogue number, lot number and expiration date.
Sterilization	Gamma irradiation between 25 kGy and 40 kGy
Component Material Toxicity	All wetted components comply with USP <88> Biological Reactivity Tests for Class VI plastics
Endotoxin level	< 2.15 EU/device for all wetted components
Assembly	This product is manufactured in a clean room environment complying with class 8 according to ISO® 14644-1.
Packaging	Each Single use pre-loaded holder is double packed.

Ordering Information

Sampling Type	Sample Volume (mL)	Sample Unit/ Single Use Holder	Qty/pack	Catalogue No.
NovaSeptum® AV Accurate Volume Sampling Unit	20mL	5 Single Sampling Unit	4	5SU461-90020
	5mL	5 Single Sampling Unit	6	5SU461-90005
NovaSeptum® High Purity Sampling Unit	50mL	5 Single Sampling Unit	8	5SU711-10050
	100mL	5 Single Sampling Unit	8	5SU711-10100
	250mL	5 Single Sampling Unit	7	5SU711-10250
NovaSeptum® Bottle Sampling Unit	1000mL	5 Single Sampling Unit	6	5SU711-11000
	60mL	5 Single Sampling Unit	6	5SU871-80060
	125mL	5 Single Sampling Unit	4	5SU871-80125
	250mL	5 Single Sampling Unit	3	5SU871-80250
	500mL	5 Single Sampling Unit	2	5SU871-80500

Sampling Type	Sample Volume (mL)	Sample Unit/ Single Use Holder	Qty/pack	Catalogue No.
NovaSeptum® AV Accurate Volume Sampling Unit	20mL	9 Single Sampling Unit	2	9SU461-90020
	5mL	9 Single Sampling Unit	3	9SU461-90005
	50mL	9 Single Sampling Unit	3	9SU711-10050
NovaSeptum® High Purity Sampling Unit	100mL	9 Single Sampling Unit	3	9SU711-10100
	250mL	9 Single Sampling Unit	3	9SU711-10250
	1000mL	9 Single Sampling Unit	3	9SU711-11000
NovaSeptum® Bottle Sampling Unit	60mL	9 Single Sampling Unit	3	9SU871-80060
	125mL	9 Single Sampling Unit	2	9SU871-80125
	250mL	9 Single Sampling Unit	1	9SU871-80250

To place and order or receive technical assistance

Find contact information for your country at:
www.emdmillipore.com/offices

For Technical Service, please visit:
www.emdmillipore.com/techservice

Merck KGaA
 Frankfurter Straße 250
 64293 Darmstadt, Germany
merckgroup.com/life-science

