

If they are precious to you, test them for *Campylobacter*

The Singlepath® Direct Campy Poultry Kit
for easy and fast on-site detection in poultry



Don't chicken out of *Campylobacter* testing

No pathogen causes human disease as often as *Campylobacter*. It can spread to an entire flock in a few days, making it a gamble to rely on the results of earlier screening tests.

Thanks to innovative Singlepath® Direct Campy Poultry Kit from Merck Millipore you can now screen shortly before slaughter because it gives you the results in less than two hours. Separating *Campylobacter* high (>10⁷ CFU per g cecal material) and low risk flocks for slaughter is a financially important decision – take it on the basis of up-to-date information!



Interested in our new test?
www.merckmillipore.com/poultry



easy

It's so simple – anyone can perform the test. All it needs is a full scoop of fresh feces. No prior enrichment step, no need for expensive lab instruments: The kit contains everything you need for testing.



fast

From sample collection to test result in less than 2 hours. The “pregnancy test” – like device displays a clear yes/no result and includes a control reaction (C) for extra reliability.



positive



negative



on-site testing

Convenient do-it-yourself test for use anywhere – even on the farm! No need to wait days for the results after you have shipped the sample to a laboratory for a real-time PCR test or for cultivation – in fact, inadequate transportation conditions can harm the sample and lead to a false result.

Ordering information

Product	Content	Ord. No.
Singlepath Direct Campy Poultry Kit	10 x Singlepath® Test Devices 10 x Dilution Tubes, 20 ml 10 x Sample Tubes, 1.5 ml 10 x Transfer Pipettes, 3 ml 1 x 1 ml Sample Buffer 1 x 100 ml Sample Diluent	1.04155.0001

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

For more information:
www.merckmillipore.com/poultry

Find contact information for your country at:
www.merckmillipore.com/offices



Merck KGaA
Frankfurter Straße 250
64293 Darmstadt, Germany
e-mail: mibio@merckgroup.com
www.merckmillipore.com/biomonitoring

Merck Millipore and the M logo are trademarks of Merck KGaA, Darmstadt, Germany.
Singlepath is a registered trademark of Merck KGaA, Darmstadt, Germany.
All other trademarks are the property of their respective owners.
Lit No. PB4377EN00 09/2013
© 2013 Merck KGaA, Darmstadt, Germany. All rights reserved.