



MOUSE ANTI-HUMAN NEUTROPHIL ELASTASE MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB1056
LOT NUMBER:	
QUANTITY:	100 µL
SPECIFICITY:	Reacts specifically with human neutrophil elastase and not with the 3 homologous neutrophil neutral serine proteases, cathepsin G, esterase N and proteinase 3. Inhibits the elastinolytic activity of elastase. Elastase is felt to play an important role in a number of disease states.
CLONE NAME:	AHN-10
ISOTYPE:	IgG ₁
FORMAT:	Ascites fluid.
STORAGE/HANDLING:	Maintain in convenient aliquots at -20°C or colder for up to 12 months. Avoid thawing and refreezing.
REFERENCES:	1. Skubitz, K.M., et al., J. Leuk. Biol. (1988), 44 :158-165.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com