



MOUSE ANTI-MEASLES (MATRIX PROTEIN) MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB8910
LOT NUMBER	
QUANTITY:	100 µg
CONCENTRATION:	1 mg/mL
SPECIFICITY:	Recognizes the measles matrix protein by indirect immunofluorescence.
APPLICATIONS:	Indirect Immunofluorescence Final working dilutions must be determined by end user.
FORMAT:	Purified immunoglobulin in 0.02M PBS with 0.25M NaCl containing 0.1% sodium azide, pH 7.6.
STORAGE/HANDLING:	Maintain at 2-8°C in aliquots for up to 12 months.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.