



HUMAN C-REACTIVE PROTEIN (CRP) PURIFIED PROTEIN

CATALOG NUMBER:	AG723
LOT NUMBER:	
QUANTITY:	1 mg
CONCENTRATION:	xx mg/mL
SOURCE:	Human pleural fluid, negative for HBsAg, anti-HCV and anti-HIV-1/2 by FDA approved assays.
PURITY:	>98% by SDS-PAGE. Single band at 26 kDa.
IEP:	No precipitin bands when tested at 2 mg/mL against anti-whole human serum.
FORMAT:	Purified protein.
PRESENTATION:	Liquid in 100mM NaCl, 10mM Tris-HCl with 0.1% sodium azide, 2mM Ca ²⁺ .
STORAGE/HANDLING:	Maintain under sterile conditions at 2-8°C for up to 1 year after date of receipt. <u>Do not freeze.</u> Full safety precautions should be taken as in the handling of any potentially infective body fluid.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.