

**MOUSE ANTI-HUMAN INTEGRIN  $\alpha$ -V  
MONOCLONAL ANTIBODY**

**CATALOG NUMBER:** MAB1978

**LOT NUMBER:**

**QUANTITY:** 100  $\mu$ L

**SPECIFICITY:** Reactive with integrin  $\alpha$ -V.

**ISOTYPE:** IgG<sub>1</sub>

**CLONE:** LM142

**APPLICATIONS:** Immunohistochemistry (Acetone Fixed Only) 1:1000  
Immunofluorescence and immunoprecipitation.  
MAB1978 is non-blocking.  
Final working dilutions must be determined by end user.

**FORMAT:** Ascites fluid with sodium azide

**STORAGE/HANDLING:** Maintain at -20°C for up to 12 months. Avoid repeated freeze/thaw cycles.

**REFERENCES:** *PNAS. USA* **84**:6471-6475 (1987).  
*Mol. Hu. Repro.* **2(7)**: 527-534 (1996).

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200  $\mu$ l or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

*For research use only; not for use as a diagnostic.*

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233  
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide  
[www.chemicon.com](http://www.chemicon.com) • [custserv@chemicon.com](mailto:custserv@chemicon.com) • [techserv@chemicon.com](mailto:techserv@chemicon.com)