

XXXXXX

7.76037.0005 01/2025

Supelco®

sigmaaldrich.com

50 Tests

1.16977

Reflectoquant®

Ammonium Test for RQflex® 20


Ammonium-Test für RQflex® 20

20 - 180 mg/l NH_4^+
Reaction time: 15 sec



Warning

Store cool and dry • (15 - 25 °C) • Kühl und trocken lagern



Please consult Instruction for Use for further information and observe the RQflex® manual.

For IfU please scan code with My M Safety App or go to sigmaaldrich.com



©2025

Lot

min. shelf life (YYYY/MM/DD)

Merk KGaA, 64271 Darmstadt, Germany
EMD Millipore Corporation, Burlington MA 01803, USA
MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville, Ontario, L6H 6J8, Canada