

## Quality Control Ranges

MILLIPLEX® Human Complement C5b-9

Catalog Number HC5B9-140K

QC Catalog Number/ Lot Number HC5B9-6140/ HC5B9-103 and HC5B9-203

<b><u>Analyte</u></b>	<b><u>QC Level</u></b>	<b><u>Expected Range</u></b>	<b><u>Units</u></b>
C5b-9	Control 1	6.3 – 13	ng/mL
	Control 2	24 – 51	ng/mL

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

The vibrant M, Milliplex, Millipore and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.

© 2026 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Cat#: HC5B9-6140, Lot#: HC5B9-103.203 (Rev 09 APR 2026)

