



GOAT ANTI-INFLUENZA B VIRUS

CATALOG NUMBER:	AB1058
LOT NUMBER:	XXXX
QUANTITY:	1 mL
CONCENTRATION:	4-5 mg/mL
PURITY:	>95%
SPECIFICITY:	Influenza B, purified virions. Non reactive with monkey kidney cells.
IMMUNOGEN:	Yamagata strain of Influenza B.
APPLICATIONS:	Applications include ELISA, fluorescence microscopy, immunoblotting; unboiled samples in 1mM DTT and 2% SDS only, immunohistochemistry. IFA: >1:1,000.
CROSS REACTIVITIES:	No cross-reactivity with Influenza A virus, Parainfluenza virus (types 1-3), RSV or Adenovirus.
FORMAT:	IgG fraction
PRESENTATION:	Liquid in PBS (0.01 M, pH 7.2) and 0.1% sodium azide.
STORAGE/HANDLING:	Maintain at +2-8°C for 3 months or at -20°C in aliquots for up to 12 months. Avoid repeated freeze/thaw cycles. Store under subdued light.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.