



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP TGF- β 1,2,3 Plex Magnetic Bead Panel Kit Catalog # TGFBMAG-64K-03

Control Catalog # TGFB-6064-3 Lot# TGFB-116 and TGFB-216

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
TGF β 1	Control I	153 – 318	pg/mL
	Control II	503 – 1045	pg/mL
TGF β 2	Control I	138 – 286	pg/mL
	Control II	501 – 1041	pg/mL
TGF β 3	Control I	129 – 268	pg/mL
	Control II	518 – 1075	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

TGFB-6064-3; TGFB-116.216 (19 December 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA