



## MOUSE ANTI-INFLUENZA A MONOCLONAL ANTIBODY BLEND

---

|                            |                                                                                                                                                                                       |                      |                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|-------------------------------------|
| <b>CATALOG NUMBER:</b>     | MAB8251                                                                                                                                                                               | <b>QUANTITY:</b>     | 100 µL                              |
| <b>LOT NUMBER:</b>         |                                                                                                                                                                                       | <b>EPITOPE:</b>      | Nucleoprotein                       |
| <b>CLONE NAME:</b>         | A1 & A3                                                                                                                                                                               | <b>HOST/ISOTYPE:</b> | Ms IgG2a (A1);IgG <sub>1</sub> (A3) |
| <b>SPECIFICITY:</b>        | All subtypes of Influenza A (H1N1, H2N2 and H3N2 and H5N1). Specific for the nucleoprotein antigen. No cross-reactivity seen to Influenza B, Parainfluenza 1, 2 or Adenovirus or RSV. |                      |                                     |
| <b>IMMUNOGEN:</b>          | Influenza A blend.                                                                                                                                                                    |                      |                                     |
| <b>APPLICATIONS:</b>       | Indirect immunofluorescence at 1:100-1:1000.<br>Optimal working dilutions must be determined by end user.                                                                             |                      |                                     |
| <b>SPECIES REACTIVITY:</b> | Reacts with Human. Reactivity with other species has not been confirmed.                                                                                                              |                      |                                     |
| <b>FORMAT:</b>             | Liquid                                                                                                                                                                                |                      |                                     |
| <b>PRESENTATION:</b>       | Ascites fluid containing no preservatives.                                                                                                                                            |                      |                                     |
| <b>STORAGE/HANDLING:</b>   | Maintain at -20°C in undiluted aliquots for up to 12 months after date of receipt.<br>Avoid repeated freeze/thaw cycles.                                                              |                      |                                     |

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC  
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.