



QUALITY CONTROL RANGES

MILLIPLEX® MAP Mouse Acute Phase Magnetic Bead Panel 2 Catalog # MAP2MAG-76K

Standard Catalog # MAP2-8076-2 Lot# 8762-8K

Analytes	Standard Value
Adipsin	3.39 ng/mL
AGP	0.525 ng/mL
A2M	7312 ng/mL

Analytes	Standard Value
CRP	20.3 ng/mL
Haptoglobin	292 ng/mL
SAP	4192 ng/mL

Control Catalog # MAP2-6076-2

<u>Apolipoproteins</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Adipsin	Control I	MAP2-108	0.031 – 0.065	ng/mL
	Control II	MAP2-208	0.19 – 0.40	ng/mL
AGP	Control I	MAP2-108	0.005 – 0.010	ng/mL
	Control II	MAP2-208	0.032 – 0.066	ng/mL
A2M	Control I	MAP2-108	72 – 149	ng/mL
	Control II	MAP2-208	446 – 926	ng/mL
CRP	Control I	MAP2-108	0.18 – 0.37	ng/mL
	Control II	MAP2-208	1.13 – 2.3	ng/mL
Haptoglobin	Control I	MAP2-108	2.8 – 5.8	ng/mL
	Control II	MAP2-208	19 – 39	ng/mL
SAP	Control I	MAP2-108	36 – 76	ng/mL
	Control II	MAP2-208	229 – 476	ng/mL

MAP2-6076-2: MAP2-108.208 (29 January 2025)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA