



**SHEEP ANTI-GLIAL DERIVED NEUROTROPHIC FACTOR (GDNF)
AFFINITY PURIFIED
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB5252P
LOT NUMBER:	
QUANTITY:	50 µg
SPECIFICITY:	Specific for GDNF. No cross reactivity with Neurturin.
IMMUNOGEN:	Recombinant human GDNF.
APPLICATIONS:	Western blot: 1 µg/mL Immunohistochemistry: 1 µg/mL Optimal working dilutions must be determined by end user.
SPECIES REACTIVITIES:	Human, rat and mouse. Other species have not been tested
FORMAT:	Affinity Purified Immunoglobulin
PRESENTATION:	Lyophilized. Reconstitute with 50 µL of sterile distilled water. Centrifuge to remove any residue. Contains no preservative. Glycerol (1:1) can be added for additional stability.
STORAGE/HANDLING:	Maintain lyophilized material at -20°C to -70°C (dry) for up to twelve months after date of receipt. After reconstitution maintain at -20°C in undiluted aliquots for up to six months. Avoid repeated freeze/thaw cycles.
RELATED REFERENCE:	Lin, L.F., <i>et al.</i> , <i>Science</i> (1993) 260 :1130-1132.

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.