



RABBIT ANTI-RAT CYTOCHROME P450 ENZYME CYP3A1 POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB1253
LOT NUMBER:	
QUANTITY:	100 µL
SPECIFICITY:	Reacts with rat cytochrome P450 enzyme CYP3A1 in hepatic microsomal fraction. No cross-reactivity with other cytochrome P450 enzymes including CYP3A2.
IMMUNOGEN:	Synthetic peptide
APPLICATIONS:	Western blot Immunohistochemistry on formaldehyde treated sections. Optimal working dilutions must be determined by end user.
FORMAT:	Rabbit antiserum
PRESENTATION:	Liquid. Contains no preservative.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months after date of receipt. Avoid repeated freeze/thaw cycles.
REFERENCE:	<i>Biochemical Pharmacology</i> (1995) 50 :2047-2056.
RELATED REFERENCES:	<i>Federation Proc.</i> (1987) 46 :2563-2566. <i>Adv. Enzymol.</i> (1987) 60 :35-87. <i>DNA and Cell Biology</i> (1993) 12 :1-51.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com