



**MOUSE ANTI-CHONDROITIN SULFATE PROTEOGLYCAN (CSPG)
PROTEIN CORE EPITOPE
MONOCLONAL ANTIBODY**

CATALOG NUMBER:	MAB1581
LOT NUMBER:	
QUANTITY:	100 µL
SPECIFICITY:	Chondroitin sulfate proteoglycan (CSPG), protein core epitope
IMMUNOGEN:	Feline brain proteoglycans
ISOTYPE:	IgM
CLONE:	Cat-315
APPLICATIONS:	Immunohistochemistry: 1:2,000 on 4% paraformaldehyde fixed frozen tissue. Immunocytochemistry: 1:500 on primary cultures of neurons. Immunoblot: 1:5,000 recognizes a band at 680 kDa. Immunoprecipitation Optimal working dilutions must be determined by the end user.
SPECIES REACTIVITIES:	Rat and feline.
FORMAT:	Ascites
PRESENTATION:	Liquid with 0.1% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months. Avoid repeated freeze/thaw cycles
REFERENCE:	Lander, C. et al., <i>J. Neuroscience</i> (1998) 18 :174-183.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2010: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com