



MOUSE ANTI-TAU, PHOSPHO THREONINE 231 MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB5450		
LOT NUMBER:			
QUANTITY:	100 µL		
SPECIFICITY:	Reacts with human Tau phosphorylated at threonine 231 and fetal tau. The antibody also reacts with dephosphorylated neurofibrillary tangles. MAB5450 is reactive with the Thr231 phosphorylated and diphosphorylated peptides. No reactivity with normal adult tau or with unphosphorylated or serine 235 phosphorylated protein.		
IMMUNOGEN:	Paired helical filaments tau preparation from human brain.		
ISOTYPE:	IgG ₁	CLONE:	PHF-6
APPLICATIONS:	Western blot Optimal working dilutions must be determined by end user.		
SPECIES REACTIVITY:	Human. Other species have not yet been tested.		
FORMAT:	Ascites fluid.		
PRESENTATION:	Liquid.		
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months after date of receipt. Avoid repeated freeze/thaw cycles.		
REFERENCE:	Bramblett, G.T., et al., <i>Neuron</i> (1993) 10 :1089-1099.		

For research use only; not for use as a diagnostic.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

For research use only; not for use as a diagnostic.

©2002 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (909) 676-9209
www.millipore.com