

don't use printout for evaluation!!

1.09535.0001 **Supelco**

**MQuant®** 100 Tests  
pH-indicator strips (non-bleeding)  
**pH 0 - 14**  
**Universal indicator**

pH-Indikatorstäbchen (nicht blutend)  
Bandelettes indicatrices de pH  
(ne déteignent pas)  
Tiras indicadoras del pH (no destiñen)

0	1	2	3	4	5	6	7
Orange	Orange	Orange	Orange	Orange	Orange	Orange	Orange
Orange	Orange	Orange	Orange	Orange	Orange	Orange	Orange
Yellow	Yellow	Yellow	Yellow	Yellow	Yellow	Green	Green
Purple	Purple	Red	Red	Red	Yellow	Yellow	Yellow

7.69735.9622 Version 07/2022 - XXXXXX

Yellow	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow	Yellow
Blue	Blue	Blue	Blue	Blue	Green	Green	Green
Red	Red	Red	Red	Red	Orange	Orange	Orange
Purple	Purple	Red	Red	Yellow	Yellow	Yellow	Yellow

14 13 12 11 10 9 8 7

Merk KGaA, 64271 Darmstadt, Germany  
Tel. +49(0)6151 72-2440  
EMD Millipore Corporation, 400 Summit Drive  
Burlington MA 01803, USA, Tel. +1-978-715-4321  
Sigma-Aldrich Canada Co. or Millipore (Canada) Ltd.  
2149 Winston Park, Dr. Oakville, Ontario, L6H 6J8  
Phone: +1 800-565-1400

Made in Germany

**pH 0 - 14**  
min. 2 s