



Quality Control Ranges
Bovine Acute Phase Protein Magnetic Bead Panel 1
Kit Cat# BAPP1-150K
Control Cat# BAPP1-6150
Lot# BAPP1-102 and BAPP1-202

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
Haptoglobin	Control 1	2.7 – 5.7	ng/mL
	Control 2	13 – 26	ng/mL
Fibrinogen	Control 1	14 – 29	ng/mL
	Control 2	64 – 132	ng/mL

BAPP1-6150: BAPP1-102.202 (27August 2025)

SigmaAldrich.com/support • SigmaAldrich.com

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.