

QUALITY CONTROL RANGES



MILLIPLEX[®] MAP Human Amyloid Beta and Tau Magnetic Bead Panel
Kit Catalog # HNABTMAG-68K
Control Catalog # HABT-6068
Lot# HABT-106 and HABT-206

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
A β 40	Control I	311 - 646	pg/mL
	Control II	1268 - 2633	pg/mL
A β 42	Control I	45 - 94	pg/mL
	Control II	181 - 375	pg/mL
tTau	Control I	175 - 364	pg/mL
	Control II	719 - 1492	pg/mL
pTau 181	Control I	12 - 24	pg/mL
	Control II	47 - 99	pg/mL

**Note: The Quality Control Ranges are generated with overnight assay format.
Quality Control values in culture media are not tested.**

HABT-6068: HABT-106.206 (23 November 2020)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.