

Don't use printout for evaluation!!

XXXXXX

7.69200.1115 11/2025

Supelco

sigmaaldrich.com

100 Tests / 25 Tests

1.10011

MQuant®

Peroxide Test Peroxid-Test



For IfU please scan code with
My M Safety App or go to
sigmaaldrich.com



0



0.5



2



5




10



25

mg/l H₂O₂



1 s



15 s

**For determination in
organic solvents see
Instruction for Use**

**Bestimmung in orga-
nischen Lösungsmitteln
s. Gebrauchsanweisung**



©2025

Store cool and dry • (15 - 25 °C) • Kühl und trocken lagern

Lot

min. shelf life (YYYY/MM/DD)

Merck KGaA, 64271 Darmstadt, Germany
EMD Millipore Corporation, Burlington, MA 01803, USA
MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville,
Ontario, L6H 5J8, Canada