



**GOAT ANTI-HUMAN IgM, AFFINITY PURIFIED
(Heavy Chain)
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AP114
LOT NUMBER:	
HOST SPECIES:	Goat
QUANTITY:	2 mg
PURITY:	Affinity Purified Antibody shows a single band against Rabbit anti Goat serum when run as the antigen in IEP at a minimum concentration of 10mg/ml. Note: The IEP also shows a second band of unknown identity.
SPECIFICITY:	Reacts with Human IgM chain. Monospecific on IEP against whole human serum. No reactivity observed against IgG and the light chains of other immunoglobulins.
CONCENTRATION:	2.0 mg/mL
FORMAT:	Nonsterile liquid. Buffer = 0.02M Phosphate Buffer, 0.25M NaCl, pH: 7.6 and 0.1% sodium azide. Filtered through a 0.45 micron membrane filter.
IMPORTANT NOTE:	Due to the wide range of variable immunoassays for which this product is intended, it is necessary that each laboratory establish individual performance characteristics for all assay components in any specific assay system.
STORAGE/HANDLING:	Store at 2-8°C.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.