

# Specification sheet

**Manganese(II) sulfate monohydrate, AR,  $\geq 99\%$**

Product number: **V800250**

CAS number: **10034-96-5**

Test	Specification
Appearance (Colour)	Pink
Appearance (Form)	Crystalline Powder
Assay (Complexometry)	$\geq 99\%$
Solubility (Colour)	Colourless
Solubility (Turbidity)	Clear
Chloride	$\leq 0.001\%$
Fe	$\leq 0.001\%$
Heavy Metals (As Pb)	$\leq 0.002\%$
pH 5 % solution in water	3.5-5