



RABBIT ANTI-HISTONE H2B POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB1623	QUANTITY:	100 µL
LOT NUMBER:			
SPECIFICITY:	Reacts with histone H2B in human nuclear fractions. No cross-reactivity with other histones, H1, H2A, H3 or H4. Binding of histones from other species has not been tested, but it is likely as the amino acid sequence of histone H2B has been found to be identical in all species examined.		
IMMUNOGEN:	Synthetic peptide.		
APPLICATIONS:	Western blot Immunohistochemistry on paraformaldehyde treated sections. Optimal working dilutions must be determined by end user.		
FORMAT:	Rabbit antiserum		
PRESENTATION:	Liquid. Contains no preservative.		
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.		

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2003 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com

AB1623/02-SEP-2011/AKA/Formatting