

1.09632.0001

Supelco

MQuant®

pH-indicator strips (non-bleeding)

pH 5.2 - 7.2

Special indicator

for pH measurement in meat

100 Tests



No contamination of meat
by soluble components.
No protein error.
Do not touch the indicator
zones.



5.2

5.6

6.0

6.4

6.8

7.2



7.69735.9619

Version 04/2025 - XXXXXX

Merck KGaA, 64271 Darmstadt, Germany
Tel. +49(0)6151 72-2440
EMD Millipore Corporation, 400 Summit Drive
Burlington MA 01803, USA, Tel. +1-978-715-4321
MilliporeSigma Canada Ltd.,
2149 Winston Park Dr, Oakville, Ontario, L6H 6J8,
Canada, Phone: +1 800-565-1400

No tocar las zonas indicadoras.
No hay error de proteínas.
No se impurifica la carne por componentes solubles.
para medir el pH en carne
Indicador especial pH 5,2 - 7,2 (no destiñen)

Ne pas toucher les zones indicatrices.
Aucune erreur de protéines.
Aucune contamination de la viande par des composants
solubles.
pour la mesure du pH dans les viandes
Indicateur spécial pH 5,2 - 7,2 (ne déteignant pas)

Indikatorzonen nicht berühren.
Bestandteile. Kein Eiweißfehler.
Keine Verunreinigung des Fleisches durch lösliche
zur pH-Messung in Fleisch
Spezialindikator pH 5,2 - 7,2 (nicht blutend)





min. 2 s



Made in Germany

pH 5.2 - 7.2