

Don't use printout for evaluation!!

XXXXXX

7.69200.2509 09/2024

Supelco

sigmaaldrich.com

100 Tests

1.10025

MQuant®

Total Hardness Test Gesamthärte-Test




©2024

°d	mg/l CaCO ₃	°e	°f
<3	<55	<4	<5
>4	>70	>5	>7
>7	>125	>9	>12
>14	>250	>18	>25
>21	>375	>26	>37

Store cool and dry • (15 - 25 °C) • Kühl und trocken lagern

Merck KGaA, 64271 Darmstadt, Germany
EMD Millipore Corporation, Burlington MA 01803, USA
MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville,
Ontario, L6H 5J8, Canada



1°d ≙ 10 mg/l CaO ≙ 1.25 °e ≙ 1.78 °f
≙ 17.8 mg/l CaCO₃ ≙ 0.18 mmol/l CaCO₃ (Ca)

Lot

min. shelf life (YYYY/MM/DD)

Download ITU
with My M
Safety

