



Quality Control Ranges
MILLIPLEX[®] MAP Mouse Myokine Magnetic Bead Panel
Catalog # MMYOMAG-74K
Control Catalog # MMY-6074 Lot# MMY-106.206

Hormone	QC Level	Expected Range	Units	Hormone	QC Level	Expected Range	Units
RBDNF	Control 1	6.7 – 14	pg/mL	MSTN	Control 1	5823 – 12094	pg/mL
	Control 2	71 – 147	pg/mL		Control 2	47023 – 97662	pg/mL
EPO	Control 1	327 – 679	pg/mL	IRISN	Control 1	688 – 1428	pg/mL
	Control 2	3084 – 6406	pg/mL		Control 2	6371 – 13231	pg/mL
IL-6	Control 1	142 – 294	pg/mL	LIF	Control 1	37 – 77	pg/mL
	Control 2	1430 – 2970	pg/mL		Control 2	353 – 733	pg/mL
IL-15	Control 1	831 – 1727	pg/mL	OSM	Control 1	235 – 487	pg/mL
	Control 2	8006 – 16628	pg/mL		Control 2	1979 – 4110	pg/mL
FGF21	Control 1	69 – 144	pg/mL	OSTCRN	Control 1	35 – 73	pg/mL
	Control 2	711 – 1477	pg/mL		Control 2	333 – 691	pg/mL
FLSTN	Control 1	2599 – 5399	pg/mL	SPARC	Control 1	2026 – 4208	pg/mL
	Control 2	24553 – 50995	pg/mL		Control 2	21236 – 44105	pg/mL

MMY-6074: MMY-106.206 (03 April 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA