

Specification sheet

Sodium hydroxide, AR, pellets, $\geq 98.5\%$

Product number: **V800383**

CAS number: **1310-73-2**

Test	Specification
Appearance (Colour)	White
Appearance (Form)	Deliquescent Pellets
Solubility (Colour)	Colourless
Solubility (Turbidity)	Clear
Assay (Acidimetry)	$\geq 98.5\%$
Carbonate	$\leq 1\%$
Fe	$\leq 0.001\%$
Pb	$\leq 0.001\%$
Potassium	$\leq 0.1\%$