



QUALITY CONTROL RANGES
MILLIPLEX® MAP Human Cardiovascular Disease Magnetic Bead Panel 2
Kit Catalog# HCVD2MAG-67K
Control Catalog # HCVMG2-6067-2 Lot# HCV2-113.213

Analyte	QC Level	Expected Range	Units	Analyte	QC Level	Expected Range	Units
ADMTS13	Control 1	18 - 37	ng/mL	NGAL	Control 1	0.18 - 0.37	ng/mL
	Control 2	136 - 283	ng/mL		Control 2	1.4 - 2.8	ng/mL
D-dimer	Control 1	3.4 - 7.1	ng/mL	sVCAM-1	Control 1	1.8 - 3.8	ng/mL
	Control 2	26 - 55	ng/mL		Control 2	14 - 29	ng/mL
GDF-15	Control 1	0.009 - 0.018	ng/mL	SAA	Control 1	2.7 - 5.6	ng/mL
	Control 2	0.07 - 0.14	ng/mL		Control 2	21 - 44	ng/mL
Myoglobin	Control 1	0.36 - 0.7	ng/mL				
	Control 2	2.5 - 5.2	ng/mL				
sICAM-1	Control 1	1.3 - 2.7	ng/mL				
	Control 2	10 - 20	ng/mL				
MPO	Control 1	0.33 - 0.7	ng/mL				
	Control 2	2.7 - 5.5	ng/mL				
P-selectin	Control 1	3.1 - 6.5	ng/mL				
	Control 2	26 - 54	ng/mL				

HCVMG2-6067-2: HCV2-113.213 (25 APRIL 2025)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA