

XXXXXX

7.76012.0007 01/2025

1.16732

Supelco®


sigmaaldrich.com

50 Tests

Reflectoquant®

Nitrite Test for RQflex® 20**Nitrit-Test** für RQflex® 200.03 - 1.00 g/l NO₂⁻

Reaction time: 60 sec



Please consult Instruction for Use for further information and observe the RQflex® manual.



For IfU please scan code with My M Safety App or go to sigmaaldrich.com



Store cold and dry • (2 - 8 °C) • Kalt und trocken lagern

Lot

min. shelf life (YYYY/MM/DD)

©2025

Merck KGaA, 64271 Darmstadt, Germany
EMD Millipore Corporation, Burlington MA 01803, USA
MilliporeSigma Canada Ltd., 2149 Winston Park Dr, Oakville,
Ontario, L6H 6J8, Canada