



RABBIT IgG PURIFIED

CATALOG NUMBER:	PP64
LOT NUMBER:	
QUANTITY:	25 mg
PROTEIN CONCENTRATION:	10.03 mg/mL
DESCRIPTION:	Purified rabbit IgG.
ELECTROPHORETIC PATTERN:	Characteristic
FORMAT:	Purified rabbit IgG.
PRESENTATION:	Frozen nonsterile liquid in a 0.01 M Phosphate buffered saline, pH 7.1. Filtered through a 0.45 micron membrane filter.
STORAGE:	Maintain at -20°C in undiluted aliquots for up to 12 months. Avoid multiple freeze-thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.



upstate | CHEMICON | *Linco*

THE EXPERTISE OF UPSTATE®, CHEMICON® AND LINCO®
IS NOW A PART OF MILLIPORE