



Parteck® LUB Excipient

A smooth operation

Achieve consistent lubrication performance.

Parteck[®] LUB is a range of vegetable-origin stearates featuring excellent batch-to-batch consistency for reliable lubrication.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

SAFC®

Pharma & Biopharma Raw Material Solutions

Parteck[®] LUB Excipient

Reliability in every sense.

When it comes to tablet production there is nothing more important than a smooth production process – in particular choosing the right lubricant. Parteck[®] LUB, a product line of our most effective excipients for lubrication, demonstrates excellent batch-to-batch consistency and helps you to achieve a reliable workflow in your tablet manufacture.

Parteck® LUB products are available in three vegetable-origin variants:

- Parteck[®] LUB MST (magnesium stearate)
- Parteck[®] LUB CST (calcium stearate)
- Parteck[®] LUB STA (stearic acid)

PARTECK[®] LUB PROVIDES:



High batch-to-batch consistency



Reliable performance thanks to specified particle size and surface area

Vegetable origin and regulatory documentation for minimal qualification and registration efforts



SEM of Parteck[®] LUB MST

Excellent batch-to-batch consistency.

To demonstrate the galenical suitability of an excipient, the lubrication effect hinges on two crucial parameters: particle size and specific surface area. We set tight limits on both. The production processes are designed to allow for consistent quality with every batch, keeping your manufacturing running smoothly.

SEM of Parteck[®] LUB CST

 $5-12 m^2/a$

BFT

Specified parameters.

We certify specific details on particle size ranges and the BET surface area for all Parteck[®] LUB. These parameters are essential in developing a stable formula for the lubricating film, enabling a reliable tableting behavior.

Vegetable origin.

As a result of the BSE/TSE crisis, many pharmacopoeias call for additional manufacturing stages to inactivate pathogens in stearates. Good news: As they are made from vegetable raw materials, our Parteck® LUB keeps your regulatory work to a minimum.

The Emprove® Program.

Ensuring the compliance of your pharma and biopharma products involves the compilation of a vast amount of data, which can be time- and resourceintensive. In order to facilitate and accelerate this process, we developed our Emprove® Program. It includes 400 pharma raw and starting materials and a selection of filtration and single-use products. Each product in the portfolio is complemented with three different types of dossiers supporting you throughout the different stages of your operations: qualification, risk assessment, and process optimization – all designed to help you speed your way through the regulatory maze.

Our knowledgeable team of technical and quality experts will share best practices and help guide you to a compliant solution that meets your performance goals.

Find out more at: **EMDMillipore.com/emprove**

	Parteck [®] LUB MST	Parteck [®] LUB CST	Parteck [®] LUB STA
d10	≥ 1 µm	≥ 0.8 mm	≥ 5 µm
d50	~ 5 µm	~ 4 µm	~ 35 µm
d90	≤ 20 µm	≤ 20 µm	≤ 200 µm

4-8 m²/a

 $0.90 - 1.40 \text{ m}^2/\sigma$

Table 1: Particle size and BET of Parteck[®] LUB.



Click. Explore. Learn more.

PARTECK® PRODUCT PORTFOLIO

Excipients for oral solid dosage forms featuring unique particle properties and outstanding individual functionalities such as solubility enhancement or controlled release. For more information, visit:

EMDMillipore.com/parteck

FORMULATION PRODUCT FINDER APP

Find the right product for specific applications at:

EMDMillipore.com/formulationapp

Ordering information

Cat. No.	Product	Pack size
1.00663.9020	Parteck® LUB MST (magnesium stearate) EMPROVE® ESSENTIAL Ph Eur, BP, JP, NF, FCC	20 kg
1.00664.9020	Parteck [®] LUB CST (calcium stearate) EMPROVE [®] ESSENTIAL Ph Eur, BP, JP, NF, FCC	20 kg
1.00661.9020	Parteck [®] LUB STA (stearic acid) EMPROVE [®] ESSENTIAL Ph Eur, BP, JP, NF	20 kg

The typical technical data above serve to generally characterize the excipient. These values are not meant as specifications and they do not have binding character. The product specification is available separately, from the website: **EMDMillipore.com**

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

For additional information, please visit **EMDMillipore.com** To place an order or receive technical assistance, please visit **EMDMillipore.com/contactPS**



MilliporeSigma, the Vibrant M, Parteck, Emprove and SAFC are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. © 2019 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Lit. No. MS_FL1198EN 09/2019