

**CHICKEN ANTI-c-myc
AFFINITY PURIFIED
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB3252
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	100 µg/mL
SPECIFICITY:	c-myc
IMMUNOGEN:	Synthetic peptide EQKLISEEDL from human c-myc.
APPLICATIONS:	Western blotting (ECL): 1:500 Immunocytochemistry: 1:250 Immunoprecipitation: 1:250 Optimal working dilutions must be determined by end user.
FORMAT:	Affinity purified chicken IgY.
PRESENTATION:	Liquid in PBS with 0.02% sodium azide, pH 7.2.
STORAGE/HANDLING:	Maintain at 2-8°C in undiluted aliquots, under sterile conditions, for up to 6 months after date of receipt. Freezing the antibody can result in loss of activity of approximately 50%. Store in the dark.

Now Available - Chicken IgY Negative Control Antibody (Cat. No. PP25),
Chicken IgY Blocking Solution (Cat. No. 2180) and Chicken IgY Precipitating Agarose
(Cat. No. 2185)

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2003 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.