

LIGHT DIAGNOSTICS™

HUMAN METAPNEUMOVIRUS (hMPV) IFA REAGENT

CATALOG NUMBER:	6002RU0
INTENDED USE:	For Research Use Only; not for use as a diagnostic.
QUANTITY:	2 mL

PRESENTATION:	One dropper vial containing a solution of mouse monoclonal antibodies in phosphate-buffered solution with Tween 20, BSA as a stabilizer and 0.1% sodium azide as a preservative.
----------------------	--

APPLICATIONS:	This reagent is "Ready to Use". No further dilution is required.
WARNINGS/PRECAUTIONS:	Sodium azide, present in this reagent, can form potentially explosive metal azides with lead and copper pipes. Upon disposal, flush with large volumes of water to prevent build-up in plumbing.
STORAGE/HANDLING:	When stored at 2-8°C, reagent is stable up to the expiration date printed on the label.

SUGGESTED PROTOCOL:	
SUBSTRATE SLIDES:	<ol style="list-style-type: none"> 1. Dissociate monolayer and wash cells 3 times with PBS. 2. Prepare a heavy cell suspension and spot onto microscope slide wells. 3. Air dry. 4. Fix in cold acetone for 10 minutes at 2°-8°C. 5. Stain immediately or store desiccated at 2°-8°C.
STAIN:	<ol style="list-style-type: none"> 1. Allow reagents to equilibrate to room temperature. 2. Add sufficient "Ready to Use" reagent to cover the cell spot. 3. Incubate for 30 minutes at 37°C in a humid chamber. 4. Rinse thoroughly with PBS/Tween 20/Azide buffer (Light Diagnostic™ Catalog No. 5015 or equivalent). 5. Add sufficient secondary FITC-labeled antibody. 6. Repeat wash step 4. 7. Mount under a coverslip using an aqueous Mounting Medium pH 8.5, (Light Diagnostic™ Catalog No. 5013 or equivalent).

Important Note: Researchers may obtain additional application information from Technical Services at: (800) 437-7500