

Specification sheet

Sodium acetate trihydrate, AR, 99-102%

Product number: **V800365**

CAS number: **6131-90-4**

Test	Specification
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Titration by HClO ₄	99-102%
pH 5 % solution in water	8.5-9.0
Chloride	≤0.0005%
Fe	≤0.0005%
Heavy Metals (As Pb)	≤0.0005%
Solubility (Colour)(C=5% solution, Water)	Colourless
Solubility (Turbidity)	Clear