



QUALITY CONTROL RANGES

HIGH SENSITIVITY GLUCAGON-LIKE PEPTIDE (GLP-1) ACTIVE CHEMILUMINESCENT ELISA CATALOG # EZGLPHS-35K

STANDARD, LYOPHILIZED CATALOG # E8035-HS

Reconstitute standard with **0.5mL** distilled or deionized water.

| | |
|----------------------|--------------------------|
| STANDARD | Active GLP1 Value |
| LOT NUMBER: 835HS-8A | 115 pM |

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6035-HS

Reconstitute each vial with **0.5mL** distilled or deionized water.

| CONTROL | LOT NUMBER | EXPECTED RANGE | UNITS |
|---------|------------|----------------|-------|
| QC 1 | GLPHS-108A | 2.6 - 5.4 | pM |
| QC 2 | GLPHS-208A | 16 - 33 | pM |

E6035-HS: GLPHS-108A.208A (30 November 2021)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA