

**SHEEP ANTI- Chx10
C-TERMINAL
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB9014
LOT NUMBER:	
QUANTITY:	250 µg
CONCENTRATION:	1 mg/mL
SPECIFICITY:	Chx10
IMMUNOGEN:	Recombinant human Chx10 protein, C-terminal.
APPLICATIONS:	Immunoblotting: 0.5-1.0 µg/mL. Reacts with the 46 kD protein in rat and mouse retinal tissue lysate. Optimal working dilutions must be determined by the end user.
SPECIES REACTIVITY:	Human, rat, mouse and bovine. Other species have not been tested.
FORMAT:	Purified immunoglobulin.
PRESENTATION:	Liquid in PBS containing 0.08% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*