



QUALITY CONTROL RANGES
Milliplex® Human Adipokine Magnetic Bead Panel 1
Kit Catalog # HADK1MAG-61K
Control Catalog # HADK1-6061-1
Lot# HADK1-111 and HADK1-211

<u>Hormone</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
Adiponectin	Control I	2633 - 5469	pg/mL
	Control II	12955 - 26906	pg/mL
Lipocalin-2/NGAL	Control I	361 - 750	pg/mL
	Control II	1684 - 3498	pg/mL
Resistin	Control I	683 - 1418	pg/mL
	Control II	3415 - 7092	pg/mL
Adipsin	Control I	1399 - 2905	pg/mL
	Control II	6044 - 12553	pg/mL
Total PAI-1	Control I	1053 - 2186	pg/mL
	Control II	5035 - 10457	pg/mL

HADK1-6061-1; HADK1-111.211 (31 January 2025)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA