Corning® Plates, Vials and Boxes

Treat Samples Safely!





Explore our promotional selection from Corning®'s expansive cell culture product range: Elplasia® Spheroid Plates enable simple, scaffold-free spheroid formation at scale for culture, imaging, and screening; Cryogenic Vials with Silicone Washers provide secure 2 mL sample storage down to −196 °C with easy identification and convenient handling; and Cryogenic Boxes, Racks & Trays offer durable, reusable, glove-friendly solutions for organized cryogenic storage.

Product Number	Product Description	Discount
Corning® Elplasia® Plates		
CLS4440-5EA	Corning® 6 Well Black Clear Round Bottom Ultra-Low Attachment Micro-Cavity plate, 500 µm Diameter and 400 µm Depth, Ind Packaged, with Lid, 1/5	25%
CLS4441-5EA	Corning® 24 Well Black Clear Round Bottom Ultra-Low Attachment Micro-Cavity plate, 500 µm Diameter and 400 µm Depth, Ind Packaged, with Lid, 1/5	25%
CLS4442-5EA	Corning® 96 Well Black Clear Round Bottom Ultra-Low Attachment Micro-Cavity plate, 500 µm Diameter and 400 µm Depth, Ind Packaged, with Lid, 1/5	25%
CLS4444-5EA	Corning® 6 Well Black Clear Square Bottom Plasma Treated, Micro-Cavity plate, 200 µm Wide and 100 µm Depth, Ind Packaged, with Lid, 1/5	25%
CLS4445-5EA	Corning® 24 Well Black Clear Square Bottom Plasma Treated, Micro-Cavity plate, 200 µm Wide and 100 µm Depth, Ind Packaged, with Lid, 1/5	25%
CLS4446-5EA	Corning® 96 Well Black Clear Square Bottom Plasma Treated,Micro-Cavity plate, 200 µm Wide and 100 µm Depth, Ind Packaged, with Lid, 1/5	25%
CLS4447-5EA	Corning® 384 Well Black Clear Square Bottom Plasma Treated, Micro-Cavity plate, 200 µm Wide and 100 µm Depth, Ind Packaged, with Lid, 1/5	25%
Screening Plates		
CLS4680-16EA	Corning® 96 Well Half Area High Content Imaging Plate, w Lid, BCB, TCT, COC, sterile	25%
CLS4517-20EA	Corning® 96 Well Half Area High Content Imaging Plate, w Lid, BCB, TCT, COC, low base, sterile	25%
CLS4681-20EA	Corning® 384 Well High Content Imaging Plate, w Lid, BCB, TCT, COC, sterile	25%
CLS4518-20EA	Corning® 384 Well High Content Imaging Plate, w Lid, BCB, TCT, COC, low base, sterile	25%
CLS4560-100EA	Corning® Microplate, 1536 Well, COC, Black Clear Bottom, NT, no Lid, Sterile, Bulk	25%
CLS4561-100EA	Corning® Microplate, 1536 Well, COC, Black Clear Bottom, TC Treated, no Lid, Sterile, Bulk	25%
CLS4563-100EA	Corning® Microplate, 1536 Well, COC, Black Clear Bottom,Corning® CellBIND® Surface, no Lid, Sterile, Bulk	25%
CLS4564-100EA	Corning® Microplate, 1536 Well, COC, Black Clear Bottom, Poly-D-Lysine, no Lid, Nonsterile, Bulk	25%

The Life Science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Product Number	Product Description	Discount
Cryogenic Vials & Boxes		
CLS431119-10EA	CRYO BOX, 81, 1–2 mL, NS,BK, 5/10	25%
CLS431120-10EA	CRYO BOX 81, 4–5 mL, NS, BK, 5/10	25%
CLS431121-10EA	CRYO BOX 100, 1–2 mL, NS, BK, 5/10	25%
CLS430488-500EA	Corning® cryogenic vial, internal thread, 2.0 mL capacity, round bottom, silicone washer seal, self-standing, polypropylene, sterile, gamma-irradiated; supplied as a case of 500. [sigmaaldrich.com]	25%
CLS431386-250EA	Corning® cryogenic vial, internal thread, 2.0 mL capacity, round bottom, silicone washer seal, self-standing, polypropylene, sterile; packaged as a case of 250. [sigmaaldrich.com]	25%
CLS430489-500EA	Corning® cryogenic vial, internal thread, 2.0 mL capacity, round bottom, silicone washer seal, polypropylene; non–self-standing version (flat packaging option). [sigmaaldrich.com]	25%
CLS431420-500EA	Corning® cryogenic vial with red cap, internal thread, 2 mL, round bottom, silicone washer, self-standing, sterile. [sigmaaldrich.com]	25%
CLS431419-500EA	Corning® cryogenic vial with green cap, internal thread, 2 mL, round bottom, silicone washer, self-standing, sterile. [sigmaaldrich.com]	25%
CLS431418-500EA	Corning® cryogenic vial with blue cap, internal thread, 2 mL, round bottom, silicone washer, self-standing, sterile. [sigmaaldrich.com]	25%
CLS431417-500EA	Corning® cryogenic vial with yellow cap, internal thread, 2 mL, round bottom, silicone washer, self-standing, sterile. [ecatalog.corning.com]	25%
CLS431421-500EA	Corning® cryogenic vial with white cap, internal thread, 2 mL, round bottom, silicone washer, self-standing, sterile. [ecatalog.corning.com]	25%
CLS431416-500EA	Corning® cryogenic vial, internal thread, 2 mL, round bottom, silicone washer seal, assorted cap colors, self-standing, sterile. [ecatalog.corning.com]	25%

Use Promotion Code ILY* to take advantage of this offer.

Discount available from 1 Oct – 31 Dec

CORNING

*Promo code cannot be combined with existing promotion

To place an order or receive technical assistance:
SigmaAldrich.com/support

For local contact information:
SigmaAldrich.com/offices

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

SigmaAldrich.com

© 2026 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_FL14618EN Ver. 1.0 71873 09/2026